**TRUSETAL
VERBANDSTOFFWERK GMBH**

LUEDDE EXOPHTHALMOMETER

#53503

The base of the transparent rule is placed directly against the external orbital canthus. The ruler must be tilted slightly upward in order to obtain the highest reading at the apex of the cornea. The patient should face good light and look directly forward. The norms are the other eye, meaning both eyes should have equal findings.

RECOMMENDATION FOR STERILIZATION OF THE LUEDDE EXOPHTHALMOMETER

The recommended solution for the tonometer sterilizer is glutaraldehyde as a mix at about 2% concentration.

Manufacturer:

MetriCide® 28 High-Level Disinfectant / Sterilant. 2.5% buffered glutaraldehyde. Recommended for immersible & heat-sensitive equipment/instruments. FDA-registered. Made by Kerr/Metrex.

Features of glutaraldehyde:

- Bactericidal, Virucidal, Fungicidal, Tuberculocidal, Sporacidal
- Convenient: no dilution necessary, just pour and use
- Noncorrosive: Compatible with all non-lensed immersible & heat-sensitive instruments/equipment
- Time-tested: Glutaraldehyde-based, 20+ year history, tried-and-true formula
- Contains surfactant: increase penetration of sterilant into instruments & equipment
- Long life: use & reuse up to 28 days

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Trusetal Verbandstoffwerk GmbH