

# Visual Development/Rehab tests

All tests and materials for visual development can also be used with elderly patients suffering from perception difficulties and in visual rehabilitation as well as for Neuro-Orthoptics.



## LEA 3-D PUZZLE®

For training and testing of eye and hand co-ordination. To be used reciprocally.

Symbol-size: 5.1 cm x 5.1 cm

Format: 17.2 cm x 17.2 cm

(GL 251600)

Stock code: 52022



## LEA MAILBOX GAME®

For assessment of visual recognition of line orientations. Patients with brain lesions may show a discrepancy between the perception of a line orientation alone (mailbox slot) and the combination with hand function (putting the "letter" in the slot).

Ø ca. 21.60 cm

(GL 254700)

Stock code: 53509



## LEA RECTANGLES GAME®

For the testing of size perception in disabled children. Contains 10 rectangles in 2 colours and 5 sizes with a surface of 10.16 cm<sup>2</sup> each.

(GL 254600)

Stock code: 53505



## LEA™ Cone Adaptation test

The test does not replace a colour perception test, but allows the assessment of whether definite colour assignment is still possible at reduced luminance levels. The test is everyday-life like and can be performed in the play situation.

Contents: 15 red, blue and white plastic chips

Format: 5 cm x 5 cm

(GL252900)

Stock code: 53506



## LEA SYMBOLS® Playing Cards

Symbols of different sizes for testing of near visual acuity in pre-school children in the play situation.

Contents: 64 cards with symbols

Card size: 4 cm x 7 cm

(GL 252500)

Stock code: 53504



## LEA™ Domino Set

The domino set allows testing near-range vision in a child-oriented manner.

Contents: 3 x 16 cards

Size: 4 cm x 7 cm

(GL251500)

Stock code: 53520



Stock code: 52096

# Visual Development/Rehab tests

All tests and materials for visual development can also be used with elderly patients suffering from perception difficulties and in visual rehabilitation as well as for Neuro-Orthoptics.



Stock code: 52099

## Form-Colour Preference test

Test to examine form and colour perception.

Contents: Set of 48 test cards, 8 motifs and frames on white and neon-coloured backgrounds (yellow, green, orange, pink, red).

Size: 25 cm x 25 cm

Stock code: 52096

## Form-Colour Preference test, version II

48 test cards with coloured figures and frames on white background. Specifically used for additional differentiated testing and therapy in subjects with perceptual problems.

Stock code: 52099



Stock code: 52021

## LEA™ Fixation Stick

With LEA SYMBOLS® and fixation targets

To effortlessly practice near vision in young children

Length: 11.5 cm

(GL 486200)

Stock code: 52021



Stock code: 52011

## LEA™ Single Presentation Cards

30 symbol cards with response plate included

Response plates included)

VA from 0.10 bis 2.50 (3 m)

Format: from 10 cm x 12.7 cm

to 12 cm x 12 cm

(GL 252700)

Stock code: 52011



Stock code:  
52052

## HEIDI® Expressions Test Game

The test enables early evaluation of vision for communication. It is particularly useful with children with multiple disabilities.

Set contains 18 cards.

Format: 10.16 cm x 10.16 cm

(GL 254505)

Stock code: 52052



Stock code:  
52026

## HIDING HEIDI® – Low Contrast Test

Playful face test for testing of low contrast vision in infants.

In the contrast levels: black, 25 %, 10 %, 5 %, 2,5 % and 1.25 %

Contents: 4 cards, printed on both sides

Format: 23 cm x 23 cm

(GL 253500)

Stock code: 52026



Stock code:  
52228

**NEW!**

## Hiding Heidi® – Low Contrast Face Test

In the contrast levels: black, 25 %, 10 %, 5 %, 2,5 % and 1,25 %

Contents: 7 cards, single printed

Format: 23 cm x 23 cm

(GL 253510)

Stock code: 52228