

OPHTHALMIC PRODUCTS

EDITION 2024/2025

EYESFIRST®

ACCESSORIES

ADHESIVE FOILS

BINOCULARITY TESTS

COLOUR TESTS

EARLY INTERVENTION

EYE DRESSINGS

LOW CONTRAST TESTS

LOW VISION TESTS

OCCLUSION EYE PATCHES

OPHTHALMOSCOPES

PRISMS

REFRACTION

STEREO TESTS

SURGICAL INSTRUMENTS

VISUAL ACUITY TESTS

VISUAL REHABILITATION



CREDITS

Internet:



www.tshs.eu

Konrad-Zuse-Str. 15 District court Bielefeld 33758 Schloss Holte-Stukenbrock HRB 34919

Germany VAT ID no.: DE 811775477

Telephone: +49 5207 991688-0 **Picture credits:** Fax: +49 5207 991688-28 © Christiane Paschke © eevl – stock.adobe.com

info@eyesfirst.eu © erikreis via Getty Images Email: orders@eyesfirst.eu © Good-Lite Company

> info@tshs.eu © HEINE Optotechnik GmbH & Co. KG

> > © Kay Pictures Ltd

© Pietrasanta Pharma S.p.A.

© Shanghai MediWorks Precision Instruments Co., Ltd.

© Shanghai VisuScience Meditech Co., Ltd.

© Zobel-Fotografie

Use of images only with Online shop: www.eyesfirst.eu prior permission of the authors © 3B Scientific

© OCULUS Optikgeräte GmbH

OPHTHALMIC PRODUCTS

Product section	Table of contents	Page
Eye Patches, Eye and Wound Dressing Materials	ORTOPAD® eye patches (SOFT boys, SOFT girls, BOYS,GIRLS, funpack, white, skin) Information on the adhesive and the material of Occlusion Eye Patches ORTOLUX® Moisture Chamber Eye Shields Eye Compresses, Spool Plasters Wound Dressing Materials	6 – 10 11 –12 14 15 16
EYESFIRST®	Table of contents	Page
Adhesive foils and Prisms	Adhesive Prism foils, Adhesive Occlusion foils, Adhesive Spherical foils Prism bars, Individual prisms, Prism Vergence sticks, Prism Vergence glasses Block Prism sets	18 19 20
Vision Tests	Paediatric Vision Tests Fixation targets, LEA GRATINGS® Tests Preferential Looking Tests Tests for Visual Development and Rehab, Cognitive LEA™ Tests Neuro-Orthoptic Tests Far distance Tests (charts) LEA™ Ring binders 3m, LEA™ Core Tests 1m Tests, Low Vision Tests, Near Vision Tests, Examination Tests KAY PICTURES, Vision Tests with borders Massachusetts Illuminated Cabinets, Translucent VA charts, Low Contrast charts ETDRS charts for Illuminated Cabinets (incl. Low Contrast charts) Low Contrast Tests (plates and ring binders incl. MARS, LEA GRATINGS®	21 22 23 24 – 25 26 – 27 28 29 – 35 36 – 37 38 – 45 46 – 51 52 – 53
Diagnostics	Binocularity Tests, Occluders, OKN Instruments, Maddox Cross Trial Frames, Trial Lenses set, Refraction and accessories, Confirmation/Prism Flippers Trial Lenses (Metal Rim), Indicators, Retina Lenses Stereo Tests (incl. LANG STEREOPAD®, TNO Test, Dinosaur Test) Colour Tests Visual Field Tests	54 - 61 62 - 65 66 - 67 68 - 69 70 - 71
Accessories	Consumables Miscellaneous – Diagnostic Eye Products	73 74 – 78
HEINE Instruments	Ophthalmoscopes Ophthalmology Sets, Retinoscopes Hand-held Instruments, Chargers Handles, Rechargeable Batteries, Spare Bulbs, Accessories	80 81 82 83 – 85
Surgical Instruments	Titanium: Forceps, Scissors, Specula, Needle Holders, Choppers, Manipulators/Hooks Stainless steel: Forceps, Scissors, Specula, Needle Holders, Cannulas Surgical Container Trays Miscellaneous (Titanium/Stainless steel) Individual surgical sets TRU-PACK®	86 – 95 96 – 100 100 101 – 102 103
Product listing	In alphabetical order	104 – 106

Expertise from a single source

Eye products TRUSETAL VERBANDSTOFFWERK GMBH

EYESFIRST®

Comprehensive portfolio for ophthalmology and eye surgery

Under the brand name EYESFIRST®, we offer diagnostic products for professional use in ophtalmology, orthoptics and optometry. EYESFIRST® is also the name of our specialist online shop. In addition, surgeons will find a large selection of surgical instruments and accessories for their surgical needs at our shop.

ORTOPAD®

Premium eye patches for children – and so much more

The development and sale of occlusion eye patches form the cornerstone of our company history. Our ORTOPAD® eye patches have been specially developed for children. Our standards in terms of material and adhesive properties, compatibility and design are equally high. In addition, we support clinics, orthoptic departments and parents with a large number of various motivational articles related to occlusion therapy and have, since 2007, been closely supporting the KARIBU project, an aid project initiated by a former colleague - who is now retired - to prevent childhood blindness in Africa.

TRU-PACK®

For more efficiency in the OR

Preparing for operations is one of the most time-consuming activities in the OR. With our tailor-made surgical sets for ophthalmology, we ensure that the surgical team has the right tools at the right time. This reduces preparation times and optimises workflows. TRU-PACK® sets contain all the necessary instruments, drapes, swabs and other accessories – all sterile and professionally packed in a practical package.



Connect with us!



We speak Social Media

Social media have fundamentally changed the communication landscape. Not only the private sector, but also companies have shifted a large part of their communication to social networks. We use our social media channels on Facebook, Instagram and Youtube to keep you up to date.

Anytime and any place

Find out what we get up to, anytime, anywhere. That way, you stay up to date with new products and promotions. You will get to know our company and our products better and can contact us quickly and easily.

Every brand has a face

A great team stands behind Trusetal Verbandstoffwerk GmbH. Get to know the people who represent our brands and our products and benefit from personal contact and exchange. We deal daily with topics you are interested in and are happy to be there for you.



@Trusetal Verbandstoffwerk GmbH@ORTOPAD eye patches

(instagram

@trusetal.verbandstoffwerk@ortopad.eyepatches

Youtube

@ORTOPADeye patches



Join our community – we are looking forward to meeting you!

KARIBU project

ORTOPAD® NEWS

KARIBU PROJECT

KENYA • CONGO • NIGERIA

We have been supporting the KARIBU project in Kenya since 2007. The aid project of our retired employee Christiane Paschke focusses on the prevention of childhood blindness through early detection and timely treatment of childhood eye diseases.

After spending several years in Kenya in a teaching position as an orthoptist, Ms Paschke initiated the project together with the Kenyan ophthalmologist Dr Margaret Njuguna. She conducts annual visits to the project sites, where she runs workshops for paediatric staff and ensures that the aid is directly benefitting the children in need.

In addition to the direct assistance we have provided for Kenya over many years, we also support Ms Paschke in her efforts to improve early detection of childhood blindness diseases in other African countries.





Our commitment in Africa for the prevention of childhood blindness



Sustainability



Packaging

We design our packaging for the ORTOPAD® mix to be as sustainable as possible. We do without unnecessary material and prefer to use bamboo cellulose packaging.



Material

When manufacturing our eye patches, we decided to use bamboo fleece. As a rapidly renewable raw material, bamboo is environmentally friendly. In addition, the fleece is very soft, robust and particularly kind to the skin.



Shipping

We do not use bubble wrap or Styrofoam for shipping. Among other things, we use packaging chips that are plant-based and thus score points with a positive environmental balance. These chips are biodegradable and can be easily disposed of with organic waste.






ORTOPAD® eye patches are a safe and gentle way to cover the eye and they will remain in place even during energetic play. They are anatomically designed with age-appropriate sizes to ensure they are an excellent fit, preventing adverse light stimuli from reaching the patched eye.



upwards of age 4 years

- skin-friendly non-woven bamboo fabric ergonomic shape with nasal slits
- breathable
- thin, soft light excluding pad
- hypoallergenic adhesive
- wide selection of motifs

Illustration of the ORTOPAD® outlines in original size:



from age 2 to 4 years

up to age 2 years





Unique texture eye patches appealing to boys with a soft fleece, ergonomic shape, nasal slits, a light-excluding centre pad and a very skin-friendly adhesive.

Stock code: IN7352	Qty: 50	67 mm x 50 mm	Junior
Stock code: IN7352	Qty: 50	76 mm x 54 mm	Medium
Stock code: IN7352	Qty: 50	85 mm x 59 mm	Regular







Unique texture eye patches appealing to girls with a soft fleece, ergonomic shape, nasal slits, a light-excluding centre pad and a very skin- friendly adhesive.

Stock code: IN73425	Qty: 50	67 mm x 50 mm	Junior
Stock code: IN73426	Qty: 50	76 mm x 54 mm	Medium
Stock code: IN73427	Qty: 50	85 mm x 59 mm	Regular







Occlusion eye patches with large scale designs for boys and girls. Ergonomic shape, nasal slits and a very skin-friendly adhesive, a lightexcluding centre pad on a soft layer ensure optimal comfort.

Junior	67 mm x 50 mm	Qty: 50	Stock code: US70237
Medium	76 mm x 54 mm	Qty: 50	Stock code: US70238
Regular	85 mm x 59 mm	Qty: 50	Stock code: US70239







The large scale designs on these occlusion eye patches were created especially to appeal to girls. Five designs with a soft layer, ergonomic shape, nasal slits, a light-excluding centre pad and a very skin-friendly adhesive.

Stock code: IN73225	Qty: 50	67 mm x 50 mm	Junior	
Stock code: IN73226	Qty: 50	76 mm x 54 mm	Medium	
Stock code: IN73227	Qty: 50	85 mm x 59 mm	Regular	







The large scale designs on these occlusion eye patches were created especially to appeal to boys. Five designs with a soft layer, ergonomic shape, nasal slits, a light-excluding centre pad and a very skin-friendly adhesive.

Junior	67 mm x 50 mm	Qty: 50	Stock code: IN73325
Medium	76 mm x 54 mm	Qty: 50	Stock code: IN73326
Regular	85 mm x 59 mm	Qty: 50	Stock code: IN73327







ORTOPAD®

White occlusion eye patches for use in occlusion treatment. ORTOPAD® is made of exceptionally comfortable and breathable fleece, is remarkably soft and especially suitable for use on all sensitive skin types. With light-excluding centre pad and very skin-friendly adhesive.

Stock code: US7002	Qty: 50	67 mm x 50 mm	Junior	
Stock code: US7002	Qty: 50	76 mm x 54 mm	Medium	
Stock code: US7002	Qtv: 50	85 mm x 59 mm	Regular	







ORTOPAD Skin

This beige-coloured occlusion eye patch is extremely gentle to skin. The soft fleece offers outstanding comfort and breathability and comes with a special light-excluding centre pad and very skin-friendly adhesive.

Junior	67 mm x 50 mm	Qty: 50	Stock code: US70145
Reaular	85 mm x 59 mm	Qtv: 50	Stock code: US70146

Information on the adhesive of ORTOPAD® Occlusion Eye Patches

PROPERTIES

The adhesive dressings of the MASTER·AID® product family, to whom ORTOPAD® Eye patches belong, are produced using an adhesive without solvent and latex free. These UV-crosslinked pressure-sensitive adhesives are applied by transfer coating and ensure excellent adhesion to the substrate. For the consumer it is especially important that the product, in comparison to other adhesives, therefore does not cause skin irritation and allergies.

APPLICATION

The pressure-sensitive adhesive is activated only by slight pressure or friction during application. It is then the body temperature which keeps the adhesive properties active. The pressure-sensitive adhesives have optimum adhesive properties, can be removed again painlessly and leave no residue on the skin.

PROTECTION OF THE ENVIRONMENT

The adhesives in use with our products ensure a clean production process. The coating does not release any solvents that are harmful for the environment or the production staff. These features ensure production of unique dressings that excel in environmental protection as well.



Bamboo as carrier material



FOR THE SAKE OF THE ENVIRONMENT

Working responsibly with nature and its raw materials is very important to us, even more so as we are a supplier of single-use medical devices. We therefore focus on sustainable, resource-saving production when choosing the materials to be used. Therefore our ORTOPAD® eye patches are made from a soft, sustainable bamboo fleece.

Bamboo scores high with its positive eco-balance. The processed raw material simultaneously has very good skin compatibility. The patches are soft, yet durable and allow the skin to breathe optimally.

With a growth rate of up to one metre per day, bamboo is one of the fastest-growing and most renewable raw materials in the world. The special thing about it is that after the harvesting of the stalks the plant remains live, because the rootstock repeatedly sprouts anew. This root system with its widely branching runners also prevents soil erosion.

Compared to traditional woods, bamboo produces much more biomass, thanks to its rapid growth with comparatively low energy input. Bamboo also contributes naturally to reduction of the greenhouse effect. Thanks to its rapid growth, bamboo absorbs a considerable amount of CO₂, whilst producing large amounts of oxygen.

The versatile applicability ultimately make it an ecologically valuable alternative to a variety of resource-scarce raw materials. Therefore we have selected bamboo as the fleece for our eye patches.



The ORTOPAD selection











ORTOLUX®

Sterile, ergonomically contoured self-adhesive eye shield with transparent, unbreakable chamber (11 mm depth). Hypo-allergenic due to a very skin-friendly adhesive. Especially suitable for maintaining a moisture chamber in cases of lagophthalmos and after care for LASIK surgeries. The shield enables a wide visual field.

Chamber size: 47 mm x 30 mm (small), 56 mm x 37 mm (large)

•

 96 mm x 68 mm (small)
 Qty: 20
 Stock code: 70106

 114 mm x 80 mm (large)
 Qty: 20
 Stock code: 70108



ORTOLUX®AÎR

Sterile, ergonomically contoured self-adhesive eye shield with transparent, unbreakable chamber (11 mm depth) and with peripheral perforation for comfortable aeration of the eye area, without visual field limitation. Especially suitable for post-operative use (e.g. after cataract operations) and for protective care of the eye.

Chamber size: 47 mm x 30 mm (small), 56 mm x 37 mm (large)

96 mm x 68 mm (small)	Qty: 20	Stock code: 70136
114 mm x 80 mm (large)	Qty: 20	Stock code: 70138



Eye shields

Aerated eye shields, for protective care of the eye and after eye-surgery, transparent.

80 mm x 66 mm Qty: 125 **Stock code: 59012**



OPTOMED®

Sterile, self-adhesive eye compresses. Made of comfortable and breathable non-woven fleece with a thick wound pad and a remarkably gentle adhesive. Suitable for quick and sterile wound care in medical practices and clinics or for after care by patients at home. Also as occluder for visual field examination.

Wound cover: 66 mm x 37 mm

96 mm x 66 mm	Qty: 10	Stock code: 70118
96 mm x 66 mm	Qty: 25	Stock code: 70119



OPTOMED® SIMPLEX

Sterile, fused edge, non-adhesive eye compresses with outstanding moisture absorbency and a thick wound pad. Also comfortable when wearing glasses and for children. For use in post-operative treatment, with corneal injuries and ointment dressings.

65 mm x 57 mm Qty: 6 **Stock code: 70204**



ROLLPORE®

White, remarkably soft plaster for highly sensitive skin; with hypoallergenic adhesive. Rolled up on an x-ray translucent plastic ring. The supporting material consists of thin polyester viscose fleece.

5 m x 1.25 cm	Qty: 18	Stock code: 70055
5 m x 2.5 cm	Qty: 12	Stock code: 70056
5 m x 5 cm	Qty: 6	Stock code: 70057



ROLLSILK®

White silk plaster with hypo-allergenic adhesive rolled up on an x-ray translucent plastic ring. Breathable, easy to tear material.

5 m x 1.25 cm	Qty: 18	Stock code: 70050
5 m x 2.5 cm	Qty: 12	Stock code: 70051
5 m x 5 cm	Qty: 6	Stock code: 70052



INJECTION und INJECTION COLOR

All around adhesive injection plaster on a roll. Breathable strip made of soft non woven fleece. Particularly indicated for wound dressing and protection after injection. Hypo-allergenic adhesive, suitable even for sensitive skin of children and seniors.

White	39 mm x 18 mm	Qty: 100	Stock code: 70425
Ladybird beetle	39 mm x 18 mm	Qty: 100	Stock code: 70481
Spooky	39 mm x 18 mm	Qty: 100	Stock code: 70477



QUADRA®MED ROUND/SQUARE

with 0.5 % Chlorhexidindigluconat

The small, white adhesive wound dressings made of soft non-woven fabric are particularly kind to the skin. The antiseptic wound pad effectively supports the prevention and treatment of bacterial wound infections. Both formats are on a roll, hygienically protected in the dispenser box and easy to tear off perforations.

rund	Ø 25 mm	Qty: 150	Stock code: 70074
quadratisch	38 mm x 38 mm	Qty: 100	Stock code: 70075



CUTIFLEX® ROUND/SQUARE

Both small-format adhesive film dressings offer waterproof and visually discreet wound protection. They are individually sealed, rolled (with perforation for tearing off) and hygienically protected in the dispenser box. The polyurethane film allows moisture to pass through unhindered from inside outwards. Two base film layers allow wrinklefree application.

round	Ø 25 mm	Qty: 150	Stock code: 70007
square	38 mm x 38 mm	Qty: 100	Stock code: 70008



STERIGRAP®

Sterile wound closure strips made of non-elastic supporting fleece. For extremely effective wound closure and to treat lacerations without the need for suturing. Only 2.5 mm (3.5 mm) wide in the middle, ends widened to 8 mm (13 mm). With hypo-allergenic adhesive.


bone	32 mm x 8 mm	Qty: 10	Stock code: 70491
stripe	75 mm x 3 mm	Qty: 10	Stock code: 70492
stripe	75 mm x 6 mm	Qty: 6	Stock code: 70493

SURGICAL INSTRUMENTS



From page 86





Please find more information also in our online shop: www.eyesfirst.eu





Each prism foil with marking of base and power

Adhesive Prism foils

For use in diagnostics and treatment, especially for pre-operative prism compensation and medium-term use with changing angles. High image quality because of high quality plastic support material.

Available power: $1\triangle$, $2\triangle$, $3\triangle$, $4\triangle$, $5\triangle$, $6\triangle$, $7\triangle$, $8\triangle$, $9\triangle$, $10\triangle$, $12\triangle$, $15\triangle$, $17.5\triangle$, $20\triangle$, $25\triangle$, $30\triangle$, $35\triangle$

TRUSETAL Prism foils size approx. 67 mm x 67 mm

OPTEGO® Prism foils size approx. 75 mm x 66 mm

Please specify the desired power when placing your order.

Adhesive Occlusion foils



TRUSETAL OCCLUSION



The strength designations denote the potential visual acuity achievable through the film, relative to a baseline decimal visual acuity of 1.00.

For use with amblyopia and strabismus treatment as an alternative to occlusion eye patches (ORTOPAD®) and glassmounted occluders. Especially suitable for treatment of diplopia and for use as gradually reduced occlusion.

Available power:

Full occlusion: beige, neutral grey, blue (patterned). Light perception foils as: light perception, < 0.10, 0.10, 0.20, 0.30, 0.40, 0.60, 0.80 The grades indicate the possible visual acuity through the film (for a reference visual acuity of 1.00 decimal).

Size:

Individual foil sizes blue: approx. 4.3 cm x 5.2 cm Individual foil sizes (except blue): approx. 5.3 cm x 5.2 cm Sheet size: approx. 21.5 cm x 10.5 cm

Please specify the desired power when placing your order.



Adhesive Spherical foils segments

For medium-term use with accommodative strabismus and hypoaccommodation prior to prescription of bifocal glasses. Also suitable for use in treatment of Low Vision.

Available power: +1.25 D, +1.50 D, +1.75 D, +2.00 D, +2.50 D, +3.00 D

Size: approx. 30 mm x 17 mm Each pack contains two lens segments.

Please specify the desired power when placing your order.



Prism Bars

Prism Bar horizontal (Luneau)

Available power: $1\triangle$, $2\triangle$, $4\triangle$, $6\triangle$, $8\triangle$, $10\triangle$, $12\triangle$, $14\triangle$, $16\triangle$, $18\triangle$,

 $20\triangle$, $25\triangle$, $30\triangle$, $35\triangle$, $40\triangle$ $45\triangle$ as individual prism.

Prism Bar vertical (Luneau)

Available power: $1\triangle$, $2\triangle$, $3\triangle$, $4\triangle$, $5\triangle$, $6\triangle$, $8\triangle$, $10\triangle$, $12\triangle$, $14\triangle$,

 $16\triangle$, $18\triangle$, $20\triangle$, $25\triangle$ $30\triangle$ as individual prism.


Prism Bar set (Luneau)

Horizontal and vertical

Power from 1 \triangle to 40 \triangle (horizontal), 1 \triangle to 25 \triangle (vertical) 30 \triangle und 45 \triangle

as individual prism.

Individual Prisms



Individual prisms rounded (Luneau)

For individual choice.

Available power:

 $2\triangle$, $4\triangle$, $6\triangle$, $10\triangle$, $12\triangle$, $14\triangle$, $20\triangle$, $50\triangle$

Individual prisms square (Gulden)

For individual choice.

Available power:

 $1\triangle,2\triangle,3\triangle,4\triangle,6\triangle,8\triangle,10\triangle,12\triangle,14\triangle,15\triangle,16\triangle,18\triangle,$

 $20\triangle$, $25\triangle$, $30\triangle$, $35\triangle$, $40\triangle$, $45\triangle$, $50\triangle$


For easy vergence testing and treatment.

For easy vergence testing and treatment.

Prism Vergence sticks $2\triangle$ and $6\triangle$

Please specify the desired power when placing your order.

Prism Vergence sticks and glasses Prism Vergence sticks 3△ and 12△ Sticks 3△ and 12△



#53024

Prism Vergence glasses
Prism power: 3△ and 12△
Without handle

Stock code: 53011

Stock code: 55401

Stock code: 55402

Stock code: 55400

Stock code:

on request

Stock code:

on request

(GL 477700)

Stock code: 53027

(GL 477702)

Stock code: 53024

(GL 477701)

GL = GOOD-LITE® stock code PRISMS



Block Prism sets

Block Prism set 22, square (Luneau)

Available power:

 $1/2\triangle, 1\triangle, 2\triangle, 3\triangle, 4\triangle, 5\triangle, 6\triangle, 7\triangle, 8\triangle, 9\triangle, 10\triangle, 12\triangle, 14\triangle, 16\triangle, 18\triangle, 20\triangle, 25\triangle, 30\triangle, 35\triangle, 40\triangle, 45\triangle, 50\triangle$ separate red lens

Block Prism set 22 rounded prisms (Luneau)

Available power:

 $1/2\triangle$, $1\triangle$, $2\triangle$, $3\triangle$, $4\triangle$, $5\triangle$, $6\triangle$, $7\triangle$, $8\triangle$, $9\triangle$, $10\triangle$, $12\triangle$, $14\triangle$, $16\triangle$, $18\triangle$, $20\triangle$, $25\triangle$, $30\triangle$, $35\triangle$, $40\triangle$, $45\triangle$, $50\triangle$ separate red lens

Block Prism set 16 square prisms (Luneau)

Available power:

 $1/2\triangle$, $1\triangle$, $2\triangle$, $4\triangle$, $6\triangle$, $8\triangle$, $10\triangle$, $12\triangle$, $15\triangle$, $20\triangle$, $25\triangle$, $30\triangle$, $35\triangle$, $40\triangle$, $45\triangle$, $50\triangle$ separate red lens

Block Prism set 16 rounded prisms (Luneau)

Available power:

 $1/2\triangle,1\triangle,2\triangle,4\triangle,6\triangle,8\triangle,10\triangle,12\triangle,15\triangle,20\triangle,\\25\triangle,30\triangle,35\triangle,40\triangle,45\triangle,50\triangle$ separate red lens

Prism set 16 pieces square (Gulden)

Available power:

 $1/2\triangle, 1\triangle, 2\triangle, 4\triangle, 6\triangle, 8\triangle, 10\triangle, \\ 12\triangle, 15\triangle, 20\triangle, 25\triangle, 30\triangle, 35\triangle, 40\triangle, 45\triangle, 50\triangle \\ \text{plus 1 red lens}$

Prism set 22 pieces square (Gulden)

Available power:

 $1/2\triangle, 1\triangle, 2\triangle, 3\triangle, 4\triangle, 5\triangle, 6\triangle, 7\triangle, 8\triangle, 9\triangle, 10\triangle, 12\triangle, 14\triangle, 16\triangle, 18\triangle, 20\triangle, 25\triangle, 30\triangle, 35\triangle, 40\triangle, 45\triangle, 50\triangle$ plus 1 red lens

Stock code: 55403

Stock code: 55405

Stock code: 55404

Stock code: 55406

Stock code: 55428

Stock code: 55427

(GL 444900)

(GL 444500)

#55428

Paediatric Vision Tests

FAR DISTANCE TESTS

LEA SYMBOLS® – single bookStock code: 52008Page 29LEA SYMBOLS®/Tumbling E Wallchart 15 linesStock code: 52035Page 28LEA SYMBOLS® – Foldable chart 10 linesStock code: 52054Page 28



NEAR VISION TESTS

LEA SYMBOLS® card

Stock code: 52016

Page 33

LEA SYMBOLS® / LEATM NUMBERS® card

Stock code: 52018

Page 33



STEREO TESTS

LANG Test I-R

LANG Test II

Stock code: 54000-R

Page 68

LANG Test II

Stock code: 54001

Page 68

LANG Test II-R

Stock code: 54001-R

Page 68

Stereo test "Fly" with LEA SYMBOLS®

Stock code: 54006

Page 69



COLOUR TESTS

Ishihara colour test14 platesStock code: 56010Page 70Waggoner colour testStock code: 56015Page 71



OPAQUE GLASS OCCLUDERS

Tiger opaque glass occluder Stock code: 59537 Page 58
Flip-Up Occluder Glasses (frosted lens) Stock code: 59598 Page 58
Flip-Up Occluder Glasses (opaque) Stock code: 59599 Page 58



OPHTHALMOSCOPES

HEINE K 180° Ophthalmology set

HEINE mini 3000° LED-Ophthalmoscope

Stock code: 57567

Page 81

Page 81



PAEDIATRIC VISUAL ACUITY SET

Consisting of #52018, #52035, #56015, #59537 **Stock code: 52233 Page 35**





Please visit our online shop: **www.eyesfirst.eu**



Fixation Targets

LANG Fixation Stick reciprocally printed, Format: 14.7 cm length

LANG Fixation Cube fixation cube with 5 objects for children

LEA™ Fixation Cube with LEA™-symbols, black/white

(GL 256400)

Stock code: 59513

Stock code: 59518

Stock code: 59530

Stock code: 59531

LEA™ Fixation Cube figures, multicoloured, red handle

(GL 256300)

GOOD-LITE® Fixation Cube motifs suitable for children, red handle

Stock code: 59595 (GL 731010)

HOTV Fixation Cube black handle

Stock code: 59565

(GL 731011)

Bernell Fixation Stick with red-green suppression test

Stock code: 59579

(GL 542050)

Heidi®/fixation paddle, medium VA: 0.50 cpcm, Format: Ø 13 cm

(GL 253100)

Heidi®/fixation paddle, small

VA: 1.00 cpcm, Format: Ø 5 cm

Stock code: 52025

Stock code: 52023

(GL 253300)

Stock code: 52024

(GL 253000)



LEA GRATINGS® tests

LEA GRATINGS® - Paddles

A preferential looking test to determine detection acuity in infants, children and adults with perception difficulties. Contains 4 paddles incl. case.

VA 0.25 cpcm, 0.50 cpcm, 1.00 cpcm, 2.00 cpcm, 4.00 cpcm and 8.00 cpcm

Format: Ø 20.3 cm

LEA GRATINGS® - Contrast Test

For testing of contrast sensitivity under different contrast levels based on 4 grades (Contrast 10 %: 0.50 cpcm, 8.00 cpcm; Contrast 2.5 %: 0.50 cpcm and 4.00 cpcm). The test serves to determine the distance in which the patient can identify the direction of lines. Includes recording forms

Format: Ø 20.3 cm

Stock code: 52051 (GL 253400)

#52051

Preferential Looking tests



For the testing of visual acuity of non-verbal children or persons with reduced vision and sensory disorders. Based on the principle of preferential looking.

Teller Acuity Cards® II

For the testing of VA in non-verbal children or in subjects with reduced vision and perceptual disorders, with central monitoring hole. For testing in 1 m / 50 cm.

Contents: 16 plates, Size: 56 cm (width) x 25 cm (height)

Teller Acuity Cards® II

Half set with 8 plates For testing in 1 m / 50 cm.



Cardiff Acuity Cards – Standard Set (for age 6 months and older)

Visual acuity test based on the principle of preferential looking. Contents: 36 plates (B to M), Size: 23 cm (width) x 30 cm (height) For testing in 1 m / 50 cm.

Cardiff Acuity Cards – Screening Set

Contents: 24 plates (LVA to M), Size: 23 cm \times 30 cm For testing in 1 m / 50 cm.

Cardiff Acuity Cards – Low Vision Set

Contents: 36 plates (LVA to J), Size: 23 cm x 30 cm For testing in 1 m / 50 cm.

Cardiff Low Vision supplementing Cards

3 additional plates for Low Vision, VA 0.03, 0.039 and 0.05 For testing in 1 m / 50 cm.



Cardiff Contrast Sensitivity Test

For the evaluation of contrast sensitivity with toddlers, infants, older children and persons with intellectual and perception impairment. All plates have images of the same size but in 12 contrast levels. 3 plates each per contrast level are based on the preferential looking principle. Contrast levels from 46 % to 1 %, equivalent to a contrast sensitivity of 2.17 to 100.

Contents: 36 plates, Format: 23 cm x 30 cm For testing in 1 m / 50 cm. Stock code: 52001

Stock code: 52001-H

Stock code: 52000

Stock code: 52092

Stock code: 52093

Stock code: 59805

Stock code: 52200

Visual Development/Rehab tests



All tests and materials for visual development can also be used with adults suffering from perception difficulties and in visual rehabilitation as well as for Neuro-Orthoptics.

LEA MAILBOX GAME®

For assessment of visual recognition of line orientations. Patients with brain lesions may show a discrepancy between the perception of a line orientation alone (mailbox slot) and the combination with hand function (putting the "letter" in the slot).

Ø approx. 21.60 cm



LEA RECTANGLES GAME®

For the testing of size perception in disabled children. Contains 10 rectangles in 2 colours and 5 sizes with a surface of 10.16 cm2 each.

(GL 254600)

(GL 252900)

Stock code: 53505

Stock code: 53506

Stock code: 53504

(GL 252500)

Stock code: 53509

(GL 254700)



LEA™ Cone Adaptation test

The test does not replace a colour perception test, but allows the assessment of whether definite colour assignment is still possible at reduced luminance levels. The test simulates normal viewing conditions and can be played as a game.

Contents: 15 red, blue and white plastic chips

Format: 5 cm x 5 cm



LEA SYMBOLS® Playing Cards

Symbols of different sizes for testing of near visual acuity in preschool children in the play situation.

Contents: 64 cards with symbols

Card size: 4 cm x 7 cm

LEA™ Domino SetStock code: 53520 The domino set allows testing near-range vision in a child-oriented (GL 251500)

manner

Contents: 3 x 16 cards Size: 4 cm x 7 cm



#53520

Visual Development/Rehab tests



All tests and materials for visual development can also be used with adults suffering from perception difficulties and in visual rehabilitation as well as for Neuro-Orthoptics.

Stock code: 52022

(GL 251600)

(GL 486200)

(GL 252700)

(GL 254505)

(GL 253500)

LEA 3-D Puzzle®

For training and testing of eye and hand co-ordination. To be used

reciprocally.

Symbol size: 5.1 cm x 5.1 cm Format: 17.2 cm x 17.2 cm

LEA™ Fixation Stick Stock code: 52021

With LEA SYMBOLS® and fixation targets
To effortlessly practice near vision in young children
VA from 0.16, 0.20, 0.25, 0.33, 0.40, 0.50, 0.66, 1.00 **(40 cm)**

Length: 11.5 cm

LEA™ Single Presentation Cards Stock code: 52011

30 symbol cards with response plate included VA from 0.10 to 2.50 **(3 m)**

#52011 Format: from 10 cm x 12.7 cm to 12 cm x 12 cm

HEIDI® Expressions Test Game Stock code: 52052

The test enables early evaluation of vision for communication.

It is particularly useful with children who have multiple disabilities.

Set contains 18 cards. Format: 10.16 cm x 10.16 cm

HIDING HEIDI® – Low Contrast Test (4 cards)

Stock code: 52026

Playful face test for testing of low contrast vision in infants. In the contrast levels: black, 25 %, 10 %, 5 %, 2.5 % and 1.25 %

Contents: 4 cards, reciprocally printed

HIDING HEIDI® – Low Contrast Test (7 cards)

Stock code: 52228

In the contrast levels: black, 25 %, 10 %, 5 %, 2.5 % and 1.25 % (GL 253510)

Contents: 7 cards, single printed

Format: 23 cm x 23 cm

#52228 Format: 23 cm x 23 cm

#52026









NEURO-ORTHOPTIC TESTS

NEAR VISION TEST CHARTS

Near vision test chart 3-1 Stock code: 52095 Page 32 Near vision test chart 4-1 Stock code: 52193 Page 32

VISUAL FIELD TESTS

Damato-Campimeter Stock code: 53508 Page 72

LEA FIXATION TARGETS HEIDI®

HEIDI®/fixation paddle, format: Ø 13 cmStock code: 52024Page 22HEIDI®/fixation paddle, format: Ø 5 cmStock code: 52025Page 22

LOW CONTRAST TEST

HIDING HEIDI® Low Contrast test 4 cards Stock code: 52026 Page 25

LOW CONTRAST RING BINDERS FOR SCREENING

LEA SYMBOLS® – Low Contrast ring binderStock code: 52028Page 53LEA NUMBERS® – Low Contrast ring binderStock code: 52075Page 53SLOAN LETTERS – Low Contrast ring binderStock code: 52181Page 53

LEA GRATINGS®

LEA GRATINGS® – Paddles Stock code: 52023 Page 22

Please visit our online shop:

www.eyesfirst.eu










PREFERENTIAL LOOKING TESTS			
Teller Acuity Cards® II	Stock code: 52001	Page 23	IIII
NYSTAGMUS GLASSES			IIII
Nystagmus glasses by Frenzel	Stock code: 53516	Page 61	
REFRACTION GLASSES FOR CHILDREN			
Tiger refraction glasses +3.00 D	Stock code: 57012	Page 63	
COLOUR TEST			
LEA™ Panel 16	Stock code: 56014	Page 71	
OKN FLAGS			
OKN flag red-white	Stock code: 53502	Page 60	
SIZE PERCEPTION –			
COMMUNICATION CAPACITY			
LEA™ RECTANGLES Game®	Stock code: 53505	Page 24	
HEIDI® Expressions Test Game	Stock code: 52052	Page 25	3
ACCESSORIES			
Reversible Occluding Glasses	Stock code: 57013	Page 76	
Visual Field Fatigue Diffuser	Stock code: 53536	Page 73	

1 1000001 11000001 iococoi 2 1000000 IOCOCOL 1000001

NEW!

#52003



#52035





#52055

NEW! 0000 0000000 #52230 0000 00 00 00





00000 00000 #52227

00000

00000

Far distance tests charts

C-Test Set of 2 plates for crowding examination (OCULUS®)

Visual acuity test (5 m) containing a set of 2 plates with Landolt-rings of different optotype interspacing as well as at a short distance of 2.6" as serial optotypes (OCULUS®), VA: 0.16 to 1.25, Format: 74.5 cm x 57.6 cm

LEA SYMBOLS®/Tumbling E Wallchart 15 lines

Reciprocally printed, response plate included (symbols) VA from 0.125 to 2.50 (3 m), Format: 43 cm x 53.4 cm

LEA SYMBOLS®/Landolt-C Wallchart 15 lines

Reciprocally printed, response plate included (symbols) VA from 0.125 to 2.50 (3 m), Format: 43 cm x 53.4 cm

LEA SYMBOLS® /Tumbling E Wallchart for 4 m / 5 m

Reciprocally printed, decimal notation only VA from 0.10 to 2.00 (4 m and 5 m), Format: 57 cm x 54.5 cm

LEA SYMBOLS® /Tumbling E Foldable chart for 4 m / 5 m

Reciprocally printed, decimal notation only VA from 0.10 to 2.00 (4 m and 5 m), Format: 57 cm x 54.5 cm

LEA SYMBOLS® - Foldable chart 10 lines

VA from 0.125 to 1.25 (3 m), Format: 40.6 cm x 40.6 cm With response plate

LEA SYMBOLS® - Wall chart 10 lines

VA from 0.125 to 1.25 (3 m), Format: 40.6 cm x 40.6 cm With response plate

LEA SYMBOLS® - Foldable chart 15 lines

VA from 0.125 to 2.50 (3 m), Format: 43 cm x 53.4 cm With response plate

LEA SYMBOLS® – Foldable chart, Optotype interspacing 100 % / 25 %

Reciprocally printed, for testing of crowding-phenomenon VA from 0.10 to 2.50 (3 m), Format: 43 cm x 54 cm Response plate and cards included

LEA Numbers® – Foldable chart 15 lines

VA from 0.125 to 2.50 (3 m), Format: 43 cm x 53.4 cm

SLOAN letters – Foldable chart

VA from 0.125 to 2.50 (3 m), Format: 45.7 cm x 45.7 cm

Landolt-C – Foldable chart, Optotype interspacing 100 % / 25 %

8 positions, reciprocally printed

VA from 0.10 to 2.50 (3 m), Format: 45.7 cm x 45.7 cm

Landolt-C – Wallchart, Optotype interspacing 100 % / 25 %

8 positions, reciprocally printed

VA from 0.10 to 2.00 (3 m), Format: 47 cm x 46 cm

Please see page 52 for 1 m charts.

Stock code: 52003

Stock code: 52035

(GL 252800)

Stock code: 52045

(GL 252850)

Stock code: 52251

(GL 250225)

Stock code: 52251-A

(GL 250224)

Stock code: 52054

(GL 250250)

Stock code: 52254

(GL 250220)

Stock code: 52055

(GL 250150)

Stock code: 52230

(GL 250160)

Stock code: 52122

(GL 271100)

Stock code: 52180

(GL 735000)

Stock code: 52220

(GL 779200)

Stock code: 52227

(GL 779400)

LEATM Ring Binders 3 m

#52008 #52236

LEA SYMBOLS® – single book

With response chart and cover plate

VA from 0.10 to 2.50 **(3 m)** Format: 12.7 cm x 12.7 cm

LEA SYMBOLS® – single book Stock code: 52236

Stock code: 52008

(GL 250600)

(GL 250625)

(GL 250700)

(GL 259200)

(GL 259225)

(GL 274000)

(GL 274400)

Stock code: 52238

Contrast: 2.5 %

#52010

#52088

#52239

With response chart and cover plate

VA from 0.10 to 2.50 **(3 m)** Format: 12.7 cm x 12.7 cm

LEA SYMBOLS® crowded symbol book Stock code: 52010

With response chart without cover plate

VA from 0.32 to 1.25 **(3 m)** Format: 12.7 cm x 12.7 cm

LEA SYMBOLS® ring binder – crowding edging Stock code: 52088

With response chart without cover plate

VA from 0.20 to 1.00 **(3 m)** Format: 10.2 cm x 10.2 cm

LEA SYMBOLS® ring binder – crowding edging Stock code: 52239

Kontrast: 2.5 %

With response chart without cover plate

VA from 0.20 to 1.00 **(3 m)** Format: 10.2 cm x 10.2 cm

LEA NUMBERS® – single book Stock code: 52009

Without response chart with cover plate

VA from 0.10 to 2.50 **(3 m)** Format: 12.7 cm x 12.7 cm

LEA NUMBERS® ring binder – crowding edging Stock code: 52237

Without response chart VA from 0.20 to 1.00 **(3 m)** Format: 10.2 cm x 10.2 cm

m 0.20 to 1.00 **(3 m)**

LEATM Core Tests Basic set



191

(6)

#52009


NEW! #52237

LEA™ Core Tests – Basic set

Please ask for details. (GL 258600)

GL = GOOD-LITE®-stock code

LEATM RING BINDERS



1 m tests

LEA SYMBOLS® - Vision chart

LEA NUMBERS® - Vision chart

With cord, foldable, response plate included VA from 0.02 to 1.00 (1 m)

Format: 40.7 cm x 54.6 cm

These charts are also for use with Low Vision

Stock code: 52029

(GL 258000)

Stock code: 52030

(GL 272100)



#52048

59896

#59806

30

Low Vision tests

LEA SYMBOLS® – Low Vision book, Set (attached at the top)

LEA SYMBOLS® – Low Vision book, Set (attached at the left)

For the examination of VA for near and distance in patients with severe visual impairment. Includes occluder for children and metal frame for adults.

VA from 0.0125 to 1.25 (3 m), 0.025 to 1.25 (40 cm)

Format: approx. 26 cm x 34 cm

Stock code: 52047

(GL 527500)

Stock code: 59807

(GL 527400)

LEA NUMBERS® – Low Vision book, Set (attached at the top)

LEA NUMBERS® – Low Vision book, Set (attached at the left)

For the examination of VA for near and distance in patients with severe visual impairment. Includes occluder for children and metal frame for

VA from 0.0125 to 1.25 (3 m), 0.025 to 1.25 (40 cm)

Format: approx. 26 cm x 34 cm


Stock code: 52048

(GL 527300)

Stock code: 59806

(GL 513100)

Low Vision tests



The images show each front and back of the reading plates.

Reading plate (german text)

Reading plate (english text)

reciprocally printed (Verdana)

Both reading plates come with a distance cord

The text in all VA lines has the same features: Continuous text of

approx. 17 words, 5 single words, numbers

VA from 0.063 to 1.25 (40 cm)

Format: 38 cm x 25 cm

Amsler test ring binder (GOOD-LITE®)

7 pages: Standard (Amsler Chart), Central (Scotoma Diagonal), Color Scotoma (Red on Black), No Lines (Dots), Metamorphopsia (Parallel Lines on Black), Metamorphopsia Reading card (Added Central Lines), Juxta-central (Standard with small squares).

LEA SYMBOLS® cards

40 M/60 M, reciprocally printed

5 cards each, can be used as individual symbols or placed in a row

Size 60 M Cards: approx. 21.9 cm x 21.9 cm Size 40 M Cards: approx. 15.25 cm x 15.25 cm

LEA NUMBERS® cards

40 M/60 M, reciprocally printed

5 cards each, can be used as individual symbols or placed in a row

Size 60 M Cards: approx. 15.25 cm x 22.9 cm Size 40 M Cards: approx. 10.2 cm x 15.25 cm



#59812

#52117



#52124

The numerator (m) indicates the test distance, the denominator (M) indicates the distance at which the gap of the equivalent Landolt ring subtends 1 minute of arc. Example for Test distance 3 m: 3/40 = 0.075

Stock code: 52056 (GL 777000)

Stock code: 52057

Stock code: 59812

Stock code: 52117

Stock code: 52124 (GL 272400)


(GL 256100)


(GL 730026)


(GL 777001)

31

#52100 #52200 #52200 #52200 #52200 #52200 #52200 #52200 #52200 #52200 #52200 #52200 #52200 #52200 #52200 #52200 #52200 #52200







Near vision tests

Birkhaeuser reading plate (German text)

decimal notations from 0.10 to 1.00 (30 cm), Format: 11 cm \times 21.5 cm

Jaeger reading card (English text)

VA from J 12 to J 1 (40 cm), Format: 14 cm x 20.5 cm

Near vision test chart 3-1

Letters, LEA NUMBERS®/LEA SYMBOLS® **(40 cm)**, with distance cord. The parallel arrangement of the letters, numbers and symbols allows for easy testing of near vision acuity and determination which optotypes can be detected. The chart is particularly useful in the diagnosis of perception and communication disorders (e.g. dementia). Legend on the back.

VA from 0.05 to 2.00, Format: $43 \text{ cm} \times 19 \text{ cm}$

Near vision test chart 4-1

Letters, LEA NUMBERS®/LEA SYMBOLS® **(40 cm)**, Landolt rings, with distance cord. Reciprocally printed near vision chart with separate legend (pdf). The parallel arrangement of the letters, numbers, symbols and Landolt rings makes it easier to determine which type of optotypes can be recognised. The chart is particularly suited for diagnostic testing in the case of sensory and communication disorders (e.g. in the case of dementia).

VA from 0.05 to 2.00, Format: 33 cm x 19 cm

Tumbling E/Landolt-C near vision card

Printed on both sides, with distance cord VA from 0.05 to 2.00 **(40 cm)**, Format: $20.3 \text{ cm} \times 25.4 \text{ cm}$

SLOAN letters card

Reciprocally printed, with distance cord VA from 0.05 to 2.00 **(40 cm)**, Format: 20.3 cm x 25.4 cm

RUNGE near vision card – Letters

For rapid screening with size progression from left to right. VA from 0.05 to 1.25 **(40 cm)**, Format: 9.5 cm x 15.9 cm

RUNGE near vision card - LEA™ NUMBERS®

For rapid screening with size progression from left to right.

VA from 0.05 to 1.25 **(40 cm)**, Format: 9.5 cm x 15.9 cm

RUNGE near vision card – LEA™ SYMBOLS®

For rapid screening with size progression from left to right. VA from 0.05 to 1.25 **(40 cm)**, Format: $9.5 \text{ cm} \times 15.9 \text{ cm}$

Stock code: 52100

Stock code: 52208

(GL 119680)

Stock code: 52095

(GL 759000)

Stock code: 52193

(GL 759100)

Stock code: 52105

(GL 728000)

Stock code: 52090

(GL 729000)

. . . .

Stock code: 52106

(GL 756000)

Stock code: 52107

(GL 756100)

Stock code: 52108

(GL 756200)

Near vision tests

The VA charts #52016, #52087, #52127, #52227 and #52020 come with one side of 50 % and 25 % optotype interspace for testing of crowding phenomenon.



LEA SYMBOLS® card

Stock code: 52016

Reciprocally printed, with response plate and cord VA from 0.05 to 2.00 (40 cm), Format: 20.3 cm x 25.4 cm (GL 250800)

(GL 251000)

(GL 270900)



LEA SYMBOLS® / LEA™ NUMBERS® card

Stock code: 52018

Reciprocally printed, with response plate and cord (symbols only) VA from 0.05 to 2.00 (40 cm), Format: 20.3 cm x 25.4 cm



LEA NUMBERS® card

Stock code: 52087

Reciprocally printed, with distance cord VA from 0.05 to 2.00 **(40 cm)**, Format: 20.3 cm x 25.4 cm



HOTV crowded card

Stock code: 52020

(GL 731000) Reciprocally printed, with response plate and cord

VA from 0.05 to 2.00 (40 cm), Format: 20.3 cm x 25.4 cm



Landolt-C near card

Stock code: 52218

(GL 779000)

8 positions, reciprocally printed, with distance cord VA from 0.05 to 2.00 (40 cm), Format: 20.3 cm x 25.4 cm



SLOAN letters card

Stock code: 52198

VA from 0.05 to 1.25 **(40 cm)**, Format: 14.0 cm x 20.3 cm (GL 725000)



Near Vision Card (englisch text)

Stock code: 52223

(GL 684600)

Near Vision Card (french text)

Stock code: 52224

Combination Near Vision and Reading Card

VA from 0.20 to 1.00 **(40 cm)**, Format: 25.4 cm x 20.32 cm

(GL 658300)

(GL 274100)



#52223

#52020

LEA NUMBERS® - Near vision card

Stock code: 52127

Reciprocally printed, without response plate (40 cm and 60 cm) VA from 0.05 to 2.00, Format: 20.3 cm x 25.4 cm



Near vision tests

Reciprocally printed, without response plate

LEA SYMBOLS® - pocket card

VA from 0.05 to 2.00 (40 cm)

Format: 9.5 cm x 15.9 cm

#52017 59896 #52019

#52094

#52119

#52126

#52186

#52253

LEA NUMBERS® - pocket card Reciprocally printed VA from 0.05 to 2.00 (40 cm) Format: 9.5 cm x 15.9 cm

Reciprocally printed

Reciprocally printed

VA from 0.05 to 2.00 (40 cm) Format: 9.5 cm x 15.9 cm

Landolt-C – pocket card Reciprocally printed VA from 0.05 to 2.00 (40 cm) Format: 9.5 cm x 15.9 cm

Landolt-C – pocket card, 8 positions Reciprocally printed VA from 0.05 to 2.00 (40 cm) Format: 9.5 cm x 15.9 cm

LEA SYMBOLS®/Landolt rings – pocket card

SLOAN-letters/Landolt rings - pocket card

SLOAN-letters – pocket card, (englisch text)

Reciprocally printed, without response plate

0000 9895

Co

000

LEA NUMBERS®/Landolt rings - pocket card Reciprocally printed VA from 0.05 to 2.00 (40 cm) Format: 9.5 cm x 15.9 cm

VA from 0.05 to 2.00 (40 cm) Format: 9.5 cm x 15.9 cm

VA from 0.05 to 2.00 (40 cm) Format: 9.5 cm x 15.9 cm

#52201

Stock code: 52201 Reading chart - high/10 % contrast English (40 cm) (GL 684700) Reading chart - high/10 % contrast Stock code: 52213 (GL 658700) French (40 cm) VA from 0.05 to 2.00 (SLOAN) and 0.20 to 1.00 (text) Format: 20.4 cm x 25.5 cm

Stock code: 52017

Stock code: 52019

Stock code: 52094

Stock code: 52219

Stock code: 52119

Stock code: 52126

Stock code: 52186

Stock code: 52253

(GL 257900)

(GL 273100)

(GL 758300)

(GL 744000)

(GL 250950)

(GL 271000)

(GL 758700)

(GL 779100)



Near vision tests

C-test - Ring Binder, Haase (OCULUS®)

Particularly suitable for near vision testing in Low Vision patients. Pages with Landolt rings in wide and narrow interspacing of 2.6". Contents: 6 pages

Stock code: 52002

Stock code: 52049

Stock code: 52050

(GL 259900)

(GL 252000)



Examination Sets

LEA SYMBOLS® - Distance Screener Set

Contents: flipchart with 10 pages, response plate VA from 0.25 to 1.25 (3 m)

Format: 8.25 cm x 22.9 cm



LEA SYMBOLS® - Near Vision Screener

Contents: 1 card, reciprocally printed with response plates and cord

VA from 0.05 to 1.25 **(40 cm)** Format: 20.3 cm x 25.4 cm

HOTV Screening flip chart (Massachusetts) For Screening pre-school children from 3 years onwards Contents: 4 plates, response cards Stock code: 52204 (GL 600304)

Contents: 4 plates, response card VA from 0.50 and 0.625 **(3 m)** Format: 8.25 cm x 22.9 cm



AAPOS Vision - Screening Set

Contents: Distance visual acuity letter chart; distant visual acuity chart with LEA SYMBOLS®; black cover for visual acuity testing, alternatively paediatric glasses occluder for visual acuity testing; 3 m distance cord; answer chart, DVD (in english)

(GL 120 AAPOS)

Stock code: 52109

Stock code: 52233

Stock code: 52234



Paediatric Visual Acuity Set

Contents: LEA SYMBOLS® / LEA NUMBERS® card #52018, LEA SYMBOLS® / Tumbling E wallchart #52035, #56015, Butterfly opaque glass occluder #53009

Suitable stereo tests can be found on pages 68 and 69.



Provikit eye tests (Project Provision)

Details and composition on request

(please refer to the illustration for the possible selection)

KAY PICTURES



KAY PICTURES provides an alternative to other vision tests designed for young children. The picture charts were developed in the UK to allow easy visual acuity screening with pictures in preschool children. The revised version is suitable for screening from 18 months of age.

The test can be performed both verbally and non-verbally using the answer chart.



KAY PICTURES Single Crowded Book (3 m)

Individual pictures with LogMAR visual acuity steps from 1.00 to -0.10 and additional near test and answer chart. The pictures each have a border to test the crowding phenomenon. Altogether, the test contains six different pictures, four per visual acuity step.

Suitable for nonverbal visual acuity screening from 18 months of age.

Stock code: 52246

(KPTSINGLE)

KAY PICTURES Linear Crowded Book (3 m)

Linear test with LogMAR visual acuity steps from 0.70 to - 0.20 and additional near test and answer chart. The pictures are arranged in a row of five boxes to produce a continuous crowding effect. Four different picture arrangements are present per visual acuity step. Suitable for nonverbal visual acuity screening from 30 months of age.

Stock code: 52247

(KPTLINEAR)



KAY PICTURES Screening Book (3 m)

Compact book in DIN A5 format with individual pictures and picture series for screening. Suitable for nonverbal vision screening from 18 months of age.

LogMAR visual acuity: 0.80 to – 0.10 **(3 m)**With answer chart

Stock code: 52248

(KPTSCREEN)



#59813

@

#52248


KAY PICTURES Low Vision Book (3 m)


This book was developed specially for screening the visually impaired. It has individual picture presentations in fourteen LogMAR visual acuity steps from 1.30 to 0.00. With additional near vision test and answer chart. Suitable for nonverbal visual acuity screening from 18 months of age.

Stock code: 59813

(KPTLOW)




















Vision Tests with borders Massachusetts

LEA Symbols® / HOTV Screening Wall chart (Massachusetts)

Reciprocally printed, with response plate

(GL 250500)

VA from 10/50 to 10/12.5 (foot) (3 m)

Format: 40.7 cm x 28 cm

LEA Symbols® Screening chart 10 lines (Massachusetts) **Stock code: 52013**Reciprocally printed, with response plate

(GL 256000)

Reciprocally printed, with response plate VA from 0.10 to 1.25 **(3 m)**Format: 40.6 cm x 40.6 cm

LEA Symbole® / Screening chart 15 lines (Massachusetts

LEA Symbols® / Screening chart 15 lines (Massachusetts) **Stock code: 52014**Reciprocally printed, with response plate
(GL 257000)

VA from 0.10 to 2.50 **(3 m)** Format: 43 cm x 53.4 cm

LEA Symbols® / HOTV Screening Wall chart (Massachusetts) Stock code: 52015

(GL 257500)

Reciprocally printed, with response plate VA from 0.20 to 2.50 **(3 m)**

....

LEA Symbols® Foldable chart with 10 lines (Massachusetts) **Stock code: 52054A**Reciprocally printed, with response plate

(GL 250200)

VA from 0.10 to 1.25 **(3 m)**Format: 40.6 cm x 40.6 cm

Format: 43 cm x 53.4 cm

LEA Symbols® Foldable chart with 15 lines (Massachusetts)

Stock code: 52249

Reciprocally printed, with response plate (GL 250100)

VA from 0.10 to 2.50 **(3 m)** Format: 43 cm x 53.4 cm

HOTV - Foldable chart (Massachusetts)Stock code: 52214Reciprocally printed, with response plate(GL 600302)

Reciprocally printed, with response plate
VA from 0.10 to 2.50 **(3 m)**Format: 45.7 cm x 45.7 cm

LEA Symbols® Near Vision Screener (Massachusetts)

Stock code: 52252

Reciprocally printed, with non-verbal response (GL 258900)
VA from 0.05 to 1.25 **(40 cm)**

Format: 20.3 cm x 25.4 cm

LEA Symbols® Screening Flip Chart (Massachusetts) Stock code: 52250

Format for preschool screening at 3 m (GL 259500)
Comprises 1 practice card, 4 pages

VA from 0.50 und 0.625 **(3 m)** Format: 8.25 cm x 22.9 cm



#52197

Illuminated Cabinets

Illuminated Cabinets allow constant indirect illumination of vision tests, as recommended by the WHO.

Specifically suitable for portable use:

Illuminated Cabinet ESV 1200 (without VA charts)Stock code: 52067Cabinet with LED lights for best illumination. Compatible with all(GL ESV1200)23 cm x 35.6 cm charts. Includes storage compartment at the back

Format: 36 cm x 23.2 cm x 5.4 cm (H/W/D)

Illuminated Cabinet Quantum ESV 1500 (without VA charts) Stock code: 52197

(GL ESV1500)

Stock code: 52069

Stock code: 52104

(GL ESC2000)

(GL ESV3000)

Corresponds in dimensions to ESV 1200 light box.

Format: 36 cm x 23.2 cm x 7.8 cm (H/W/D)

Please see matching Visual Acuity charts on page 41 to page 45.

Variant to ESV 1200, calibrated for study purposes (similar to ESV 3000).



Illuminated Cabinets ETDRS: Self-calibrating ESV 3000 ETDR: Illuminated cabinet with LED light:

Self-calibrating ESV 3000 ETDRS Viewer (without VA charts) Illuminated cabinet with LED lights (can be wall mounted), illumination for both photopic and mesopic light levels.

Choice of 3 adjustable luminance levels (the viewer automatically calibrates to 85 cd/m² or 3 cd/m²).

Format: approx. 62 cm x 65 cm x 9.75 cm (H/W/D)

#52069 OF COHZV
SZNDC
VKCNR
KCNHN
ZKDVC
HVORK
HVORK
HVORK
Chart and additional glare lights

#52104



ESC 2000 ETDRS Viewer (without VA charts)

The ESC 2000 ETDRS Viewer is a competitive alternative to the ETDRS Viewer ESV 3000. All ETDRS charts fit both viewers. Format: approx. 62 cm x 65 cm x 7.65 cm (H/W/D)

Please see matching Visual Acuity charts on page 46 to 51. ETDRS sets available on request.



#52070

Illuminated Cabinets accessories

ESV 3000 Glare lights for ETDRS Viewer ESV 3000

Stock code: 52168

(GL 500100)

Floor stand ESV 3000 for Good-Lite Illuminated ETDRS cabinets

Stock code: 52070 (GL 500004)

Floor stand for ESV 1200 and ESV 1500

Stock code: 52178

Stock code: 52207

Availability

upon request

(GL 705000)





#52207



LUVIS mini - Illuminated Cabinet

LUVIS mini - Illuminated Cabinet

Standardized conditions for testing contrast vision. Appropriate for MARS contrast tests, smaller eye chart formats and colour tests with a size up to approx. 38 cm x 28 cm.

Examinations to record visual acuity and contrast threshold require homogeneous illumination of the visual tests according to DIN EN ISO 8596:2009 and DIN 58220. These requirements are not met with conventional room lighting or a desk lamp.

Luvis mini uses diffuse reflected light to allow even light distribution over the entire eye chart. An illumination difference of only 3% is made possible with the device; at present 25% tolerance in uniformity of illumination is specified.

Advantages of the LUVIS mini:

- Meets the lighting specifications of EN ISO 8596 (luminance ca. 200 cd/m²)
- Achieves standardized, comparable and reproducible results
- Innovative and simple lighting solution

Format: 47 cm (height) x 37 cm (width) x 16 cm (depth)


Eye charts sold separately




#52163



#52137



#52091



#52138



#52139

CSV System

CSV 1000 HGT viewer with glare lighting and 1 test chart #52138 The gold standard for measuring contrast vision and glare sensitivity. Self-calibrating light system.

Please see the matching CSV charts at the bottom of this page.

Weight: approx. 13.6 kg

Format: 69.69 cm x 53.34 x 10 cm

Illuminated viewer CSV 1000 and 1 test chart #52091

Standardised contrast test device with corresponding CSV vision testing chart. Self-calibrating light system without glare sensitivity testing.

Please see the matching CSV charts at the bottom of this page.

Weight: approx. 13.6 kg

Format: 69.69 cm x 53.34 x 10 cm

Stock code: 52163

(GL 500050) Availability upon request

Stock code: 52137

(GL 500000) Availability upon request

Please see page 37 for the matching floor stand.

CSV charts only suitable for CSV-1000 HGT viewer

CSV Contrast chart test "1000E"

Contrast chart with grid pattern for differentiated contrast sensitivity testing. At the recommended test distance of 2.5 m the grid patterns will test 3, 6, 12 and 18 cycles / degree.

Does not fit ETDRS Viewers. Format: $49.5 \text{ cm} \times 44.5 \text{ cm}$

CSV contrast chart "1000S"

The CSV contrast chart "1000S" is generally used for contrast screening for cataracts in the office or for mobile clinics. The uppermost segment with the usual visual acuity test series from 0.1 to 1.33, a segment with grid pattern (6 or 12 cycles/degree), lowest segment shows realistic traffic situations.

Format: $49.5 \text{ cm} \times 44.5 \text{ cm}$

CSV 1000 ETDRS chart – SLOAN letters "R"

The CSV-1000 ETDRS chart "R_{*}, offers a standardised test surface, based on American standard illumination 85 cd/m2.

Format: $49.5 \text{ cm} \times 44.5 \text{ cm}$

Stock code: 52091

(GL 500051) Availability upon request

Stock code: 52138

(GL 500001) Availability upon request

Stock code: 52139

(GL 500003) Availability upon request

C O O C O

#52058

Vision charts distance translucent

LogMAR charts for illuminated cabinets #52067 and #52197. Please note the various test distances.

Landolt-C chart, translucent

Repetitive groups

Without response plate

VA from 0.125 to 2.00 (3 m)

Format: 23 cm x 35.6 cm

Landolt-C chart, translucent Stock code: 52221

Stock code: 52058

(GL 600703)

(GL 779300)

8 positionsRepetitive groups

Without response plate

VA from 0.125 to 2.00 **(3 m)**

Format: 23 cm x 35.6 cm

LEA SYMBOLS® chart, translucent Stock code: 52059

Repetitive groups (GL 252400)

With response plate VA from 0.16 to 2.50 (3 m)

Format: 23 cm x 35.6 cm

LEA SYMBOLS®, translucent chart with 13 lines Stock code: 52112

With response plate (GL 250400)

VA from 0.16 to 2.50 **(3 m)** Format: 23 cm x 35.6 cm

LEA NUMBERS® chart, translucent Stock code: 52060

Repetitive groups (GL 271200)

Without response plate

VA from 0.16 to 2.50 **(3 m)** Format: 23 cm x 35.6 cm

Tumbling E chart, translucent Stock code: 52202

Linear-spaced Distance chart in American notation (20/20) (GL 600738)

VA from 0.20 to 1.25 **(3 m)** Format: 23 cm x 35.6 cm

Tumbling E chart, translucent wide-spaced Stock code: 52205

Linear-spaced Distance chart in American notation (20/20) and wide

inter-optotype distance


VA from 0.20 to 1.25 **(3 m)** Format: 23 cm x 35.6 cm


HOTV chart, translucent Stock code: 52062

Without response plate (GL 600736)

VA from 0.16 to 2.50 bzw. 20/125 to 20/8 **(3 m)**

Format: 23 cm x 35.6 cm






#52060

69586

68965

96859 ...



#52205



1000

MEBWE

WEBMEM

#52202

.....


#52062

(GL 600722)

#52061



#52171







#52123

Vision charts distance translucent

LogMAR charts for illuminated cabinets #52067 and #52197. Please note the various test distances.

SLOAN letters chart with 13 lines, translucent

Repetitive groups Without response plate VA from 0.16 to 2.00 **(4 m)** Format: 23 cm x 35.6 cm

SLOAN letters chart, translucent

Repetitive groups Without response plate VA from 0.12 to 2.00 **(3 m)** Format: 23 cm x 35.6 cm

SLOAN letters chart, translucent

Zufallsanordnung mit Repetitive groups Without response plate VA from 0.16 to 2.00 **(4 m)** Format: 23 cm x 35.6 cm

SLOAN letters chart with 13 lines, translucent

Repetitive groups Without response plate VA from 0.16 to 2.00 **(6 m)** Format: 23 cm x 35.6 cm

LEA SYMBOLS® - Translucent chart with 12 lines

Repetitive groups With response plate VA from 0.20 to 2.50 **(4 m)** Format: 23 cm x 35.6 cm

LEA SYMBOLS® – Translucent chart with 12 lines

Repetitive groups Without response plate VA from 0.20 to 2.50 **(4 m)** Format: 23 cm x 35.6 cm

LEA NUMBERS® – Translucent chart with 12 lines

Repetitive groups Without response plate VA from 0.20 to 2.50 **(4 m)** Format: 23 cm x 35.6 cm Stock code: 52172

(GL 600739)

Stock code: 52061

(GL 600705)

Stock code: 52240

(GL 600762)

Stock code: 52171

(GL 600737)

Stock code: 52115

(GL 255100)

Stock code: 52116

(GL 255101)

Stock code: 52123

(GL 271300)

Vision charts near translucent

LogMAR charts for illuminated cabinets #52067 and #52197. Please note the various test distances.

LEA SYMBOLS® - Translucent chart

With response plate

VA from 0.05 to 2.00 **(40 cm)** Format: 23 cm x 35.6 cm

Stock code: 52118

(GL 257800)

LEA NUMBERS® – Translucent chart near with 12 lines

Without response plate VA from 0.05 to 2.00 **(40 cm)** Format: 23 cm x 35.6 cm Stock code: 52125

(GL 273000)

#52125



SLOAN letters - Translucent near vision chart

Without response plate VA from 0.05 to 2.00 **(40 cm)** Format: 23 cm x 35.6 cm Stock code: 52184

(GL 756300)

SLOAN letters - Translucent near vision chart

Without response plate VA from 0.05 to 2.00 **(40 cm)** Format: 23 cm x 35.6 cm Stock code: 52185

(GL 756400)

#52184



#52185

#52190 #52188 #52192

Vision charts intermediate distance translucent

Stock code: 52188

Stock code: 52190

Stock code: 52192

Stock code: 52212

Stock code: 52187

Stock code: 52189 (GL 758500)

Stock code: 52191 (GL 758600)

(GL 758410)

(GL 758510)

(GL 758610)

(GL 756333)


(GL 758400)

LogMAR charts for illuminated cabinets #52067 and #52197. Please note the various test distances.

Translucent chart, near - SLOAN letters, chart "1"

Without response plate VA from 0.05 to 2.00 (66 cm)

Format: 23 cm x 35.6 cm



Without response plate VA from 0.05 to 2.00 (66 cm)

Translucent chart, near – SLOAN letters, chart "3"

Without response plate VA from 0.05 to 2.00 (66 cm)

Translucent chart, near - SLOAN letters, chart "2"

VA from 0.05 to 2.00 (80 cm)

Translucent chart, near – SLOAN letters, chart "3,,

Format: 23 cm x 35.6 cm

Format: 23 cm x 35.6 cm

Translucent chart, near – SLOAN letters

Without response plate VA from 0.05 to 2.00 (66 cm) Format: 23 cm x 35.6 cm

Translucent chart, near – SLOAN letters, chart "1"

Without response plate VA from 0.05 to 2.00 (80 cm) Format: 23 cm x 35.6 cm

Without response plate Format: 23 cm x 35.6 cm

Without response plate VA from 0.05 to 2.00 (80 cm) Format: 23 cm x 35.6 cm




#52187

ZCRDH SOK

RNOVS

#52191





#52078



#52083



#52173

Low Contrast charts translucent

LogMAR charts for illuminated cabinets #52067 and #52197. Please note the various test distances.

LEA SYMBOLS® Low Contrast chart – 1.25 % Stock code: 52076

(GL 253700)

LEA SYMBOLS® Low Contrast chart – 2.5 % Stock code: 52027

(GL 253800)

LEA SYMBOLS® Low Contrast chart – 5 % Stock code: 52077

(GL 253900)

LEA SYMBOLS® Low Contrast chart – 10 % Stock code: 52078

(GL 254000)

(GL 254100)

LEA SYMBOLS® Low Contrast chart – 25 % Stock code: 52079

Charts for **3 m** distance VA from 0.12 to 2.50 With response plate

LEA NUMBERS® Low Contrast chart – 1.25 % Stock code: 52080

(GL 271500)

LEA NUMBERS® Low Contrast chart – 2.5 % Stock code: 52081

(GL 271600)

LEA NUMBERS® Low Contrast chart – 5 % Stock code: 52082

(GL 271700)

LEA NUMBERS® Low Contrast chart – 10 % Stock code: 52083

(GL 271800)

(GL 271900)

LEA NUMBERS® Low Contrast chart – 25 % Stock code: 52084

Charts for **3 m** distance VA from 0.12 to 2.50

SLOAN letters Low Contrast chart – 1.25 % Stock code: 52177

(GL 600754)

SLOAN letters Low Contrast chart – 2.5 % Stock code: 52176

(GL 600753)

SLOAN letters Low Contrast chart – 5 % Stock code: 52175

(GL 600752)

SLOAN letters Low Contrast chart – 10 % Stock code: 52174

(GL 600751)

SLOAN letters Low Contrast chart – 25 % Stock code: 52173

Charts for **3 m** distance (GL 600750)

LOW CONTRAST CHARTS | 45

VA from 0.12 to 2.50

00000

00000

00000

00000


00000

#52063

#52167

#52216

#52215



#52217



#52159

ETDRS charts LOGMAR

ETDRS translucent charts for all ETDRS Viewers. Please note the various test distances.

ETDRS chart - Landolt-C

Translucent chart for **4 m** distance, **4 positions** VA from 0.10 to 2.00 Format: 62 cm x 65 cm

ETDRS chart - Landolt-C

Translucent chart for **4 m** distance, **8 positions** VA from 0.10 to 2.00 Format: 62 cm x 65 cm

ETDRS chart – Landolt-C, scrambled version

Translucent chart for **4 m** distance, **8 positions** VA from 0.10 to 2.00 Format: 62 cm x 65 cm

ETDRS chart – Landolt-C

Translucent chart for $\bf 3 \ m$ distance, $\bf 8 \ positions$ VA from 0.10 to 2.00 Format: 62 cm \times 65 cm

ETDRS chart – Landolt-C, scrambled version

Translucent chart for $\bf 3~m$ distance, $\bf 8~positions$ VA from 0.10 to 2.00 Format: 62 cm \times 65 cm

ETDRS chart - HOTV

Translucent chart for **4 m** distance VA from 0.10 to 2.00 Format: 62 cm x 65 cm

ETDRS chart – HOTV with 50 % Crowding

Translucent chart for **4 m** distance VA from 0.20 to 1.25 Format: 62 cm x 65 cm

ETDRS chart – Tumbling E

Translucent chart for **4 m** distance **ETDRS chart – Tumbling E**

Translucent chart for **3 m** distance VA from 0.10 to 2.00 Format: 62 cm x 65 cm

Stock code: 52063

(GL 500037)

Stock code: 52167

(GL 500070)

Stock code: 52215

(GL 500071)

Stock code: 52216

(GL 500072)

Stock code: 52217

(GL 500073)

Stock code: 52064

(GL 500025)

Stock code: 52229

(GL 500035)


Stock code: 52159

(GL 500044)

Stock code: 52243

(GL 500061)

ETDRS charts Original series LOGMAR



ETDRS translucent charts for all ETDRS Viewers. Please note the various test distances.


ETDRS Original series 4 m – SLOAN letters, chart "R"

ETDRS Original series 4 m – SLOAN letters, chart "1"

VA from 0.10 to 2.00

ETDRS Original series 4 m – SLOAN letters, chart "2" Stock code: 52098 Translucent chart for 4 m distance (GL 500015)

Format: 62 cm x 65 cm



ETDRS Original series 3 m – SLOAN letters, chart "R"

ETDRS Original series 3 m – SLOAN letters, chart "1"

ETDRS Original series 3 m – SLOAN letters, chart "2"

Translucent chart for 3 m distance

VA from 0.10 to 2.00 Format: 62 cm x 65 cm

ETDRS Original series 2 m – SLOAN letters, chart "R"

ETDRS Original series 2 m – SLOAN letters, chart "1"

ETDRS Original series 2 m – SLOAN letters, chart "2"

Translucent chart for 2 m distance

VA from 0.10 to 2.00 Format: 62 cm x 65 cm Stock code: 52160 (GL 500046)

Stock code: 52037

Stock code: 52097

(GL 500013)

(GL 500014)

Stock code: 52161

(GL 500047) Stock code: 52162

(GL 500048)

Stock code: 52153

(GL 500038)

Stock code: 52154

(GL 500039)

Stock code: 52155

Stock code: 52169

Stock code: 52170


(GL 511113)

(GL 511114)

(GL 500040)



The roll-up charts are made of thin styrene plastic material and can be easily rolled up for transport and storage. Each of these charts can be suspended with a wall hanger.



Translucent chart for **4 m** distance, can be rolled up

(letter sequence as #52160)

VA from 0.10 to 2.50, Format: 62 cm x 65 cm

ETDRS roll-up chart "1" - Original series


Translucent chart for **4 m** distance, can be rolled up

(letter sequence as #52161)

VA from 0.10 to 2.50, Format: 62 cm x 65 cm



47



ETDRS charts LEA™

ETDRS translucent charts for all ETDRS Viewers. Please note the various test distances.

ETDRS chart - LEA SYMBOLS®

Translucent chart for 4 m distance VA from 0.125 to 2.50 Format: 62 cm x 65 cm

Stock code: 52114

(GL 253333)

95685 59869 8 6 56898 6 9 65989 89 95869 96 #52065 -58 68958 96895 58659 #52195 69586-

ETDRS chart - LEA NUMBERS®

Response chart included

ETDRS chart "A" - LEA NUMBERS®

(GL 270200) Stock code: 52195 (GL 270200A)

Stock code: 52065

ETDRS chart "B" - LEA NUMBERS®

Translucent chart for 4 m distance VA from 0.10 to 2.50 Format: 62 cm x 65 cm

Stock code: 52196 (GL 270200B)

#52196

ETDRS chart – LEA NUMBERS®

Translucent chart for 3 m distance VA from 0.10 to 2.50

(GL 271400)

Stock code: 52065A

Format: 62 cm x 65 cm



Translucent chart for 3 m distance VA from 0.10 to 2.50 Format: 62 cm x 65 cm

Response chart included

Stock code: 52066

(GL 250300)

59869 86956 69585 96859-#52065A 00000 0000

95685

#52066

The roll-up charts are made of thin styrene plastic material and can be easily rolled up for transport and storage. Each of these charts can be securely suspended with a wall hanger.



ETDRS roll-up chart - LEA NUMBERS®

Translucent chart for 4 m distance

VA from 0.10 to 2.50 Format: 62 cm x 65 cm

Response chart included

Translucent chart for 3 m distance VA from 0.10 to 2.50 Format: 62 cm x 65 cm Response chart included

ETDRS roll-up chart - LEA SYMBOLS®



#52121


0000 00000 00000 00 #52111

Stock code: 52121 (GL 270211)

Stock code: 52111

(GL 250311)

ETDRS charts "2000" LogMAR



ETDRS translucent charts for all ETDRS Viewers. Please note the various test distances.

ETDRS "2000" – SLOAN letters, chart "1"

Stock code: 52140

(GL 500016)

ETDRS "2000" - SLOAN letters, chart "2"

Stock code: 52141

(GL 500017)


ETDRS "2000" – SLOAN letters, chart "3"

Stock code: 52142

Chart for 4 m distance VA from 0.10 to 2.00

Format: 62 cm x 65 cm

(GL 500018)



ETDRS "2000" - SLOAN letters, chart "1"

Stock code: 52231

(GL 500010)

ETDRS "2000" - SLOAN letters, chart "2" Stock code: 52255

(GL 500011)

ETDRS "2000" - SLOAN letters, chart "3"

Stock code: 52244

Chart for **3 m** distance VA from 0.10 to 2.00

Format: 62 cm x 65 cm

(GL 500012)



ETDRS "2000" - SLOAN letters, chart "1"

Stock code: 52156

(GL 500041)

ETDRS "2000" - SLOAN letters, chart "2"

Stock code: 52157

Stock code: 52158

(GL 500042)


(GL 500043)

ETDRS "2000" - SLOAN letters, chart "3"

Chart for 2 m distance

VA from 0.10 to 2.00

Format: 62 cm x 65 cm



#52158

49

ETDRS CHARTS

ETDRS Low Contrast charts LogMAR



#52143



#52071



#52072



#52073



#52074

ETDRS translucent charts for all ETDRS Viewers. Please note the various test distances.

ETDRS Low Contrast chart – SLOAN letters 1.25 % Stock code: 52143

STORE I and Combinate object. CLOAN letters 0.5 % Charles and a 5007

ETDRS Low Contrast chart – SLOAN letters 2.5 % Stock code: 52071

ETDRS Low Contrast chart – SLOAN letters 5 % Stock code: 52072

(GL 500020)

ETDRS Low Contrast chart – SLOAN letters 10 % Stock code: 52073

(GL 500019)

ETDRS Low Contrast chart – SLOAN letters 25 % Stock code: 52074

Charts for **4 m** distance: VA from 0.10 to 2.00 Format: 62 cm x 65 cm

ETDRS Low Contrast chart – SLOAN letters 2.5 %

ETDRS Low Contrast chart – SLOAN letters 5 %

ETDRS Low Contrast chart – SLOAN letters 10 %

ETDRS Low Contrast chart – SLOAN letters 25 %

Charts for **3 m** distance: VA from 0.10 to 2.00 Format: 62 cm x 65 cm Stock code: 52147

(GL 500029)

(GL 500023)

(GL 500022)

Stock code: 52146

(GL 500028)


Stock code: 52145

(GL 500027)

Stock code: 52144

(GL 500026)

ETDRS Low Contrast charts LogMAR



ETDRS translucent charts for all ETDRS Viewers. Please note the various test distances.

ETDRS Low Contrast chart - LEA NUMBERS® 2.5 %

Stock code: 52128

(GL 274200)

ETDRS Low Contrast chart - LEA NUMBERS® 10 %

Stock code: 52129

Chart for **4 m** distance VA from 0.10 to 2.00

Format: 62 cm x 65 cm

(GL 274300)


Evans Low Contrast charts



#52151

The Evans contrast charts are suitable for all ETDRS Viewers.

The Evans letter contrast test (ELCT) is a standardised test for testing contrast sensitivity using SLOAN letters. Displayed at a test distance of 1 m, the optotype size corresponds to a visual acuity of 0.031. Each of the 3 consecutive letters has the same contrast value. The contrast increases from a group of three to the next group of three in increments of 0.15 log.



Stock code: 52151

12 contrast levels from standard 100 % contrast to 1.65 log contrast or 2.24% contrast. Studies show that 1.65 log contrast corresponds to the contrast sensitivity of the normal elderly population (based on the Weber contrast)

(GL 500331)

Format: 62 cm x 65 cm

Evans Contrast test – ELCT chart "2"

Evans Contrast test – ELCT set

Stock code: 52152

4 contrast levels from 1.80 to 2.25 log contrast (GL 500332)

ELCT chart "1" and ELCT chart "2" as a set for testing contrast vision

Stock code: 52150 (GL 500330)

All Evans contrast charts have a format of 62 cm x 65 cm.







Low Contrast tests charts

MARS Low Contrast test – letters

MARS Low Contrast test – numbers


Three plates with 48 optotypes in 0.04 log steps Format: $22.8 \text{ cm} \times 35.6 \text{ cm}$

Stock code: 52085

Availability upon request

Stock code: 52086

Availability upon request



LEA SYMBOLS®/ E-hooks - Low Contrast Wallchart

LEA SYMBOLS®/ Landolt rings – Low Contrast Wallchart

15 lines in 2.5 % contrast, reciprocally printed, response plate included VA from 0.10 to 2.50 (3 m), Format: $43 \text{ cm} \times 53.4 \text{ cm}$

Stock code: 52036

(GL 254900)

Stock code: 52046

(GL 254950)



LEA GRATINGS® Low Contrast test

For testing of contrast sensitivity at diff erent contrast levels based on 4 grades (Contrast 10 %: 0.5 cpcm, 8.0 cpcm; Contrast 2.5 %: 0.5 cpcm and 4.0 cpcm). The test serves to determine the distance at which the patient can identify the direction of the lines. Includes recording forms

Format: Ø 20.3 cm

Stock code: 52051

(GL 253400)

Low Contrast tests ring binders





NEW!#52236 #52239



Adult Near Low Contrast Test – ring binder

7 pages with ETDRS Letters and Continuous Text in 100 %, 25 %, 10 %, 5 %, and 2.5 % contrast level.

VA from 0.20 to 1.00 as reading text, 0.05 to 2.00 as ETDRS chart **(40 cm)** Format: $15.8 \text{ cm} \times 15.8 \text{ cm}$

Paediatric Near Low Contrast Test – ring binder

7 pages in full contrast (100 %), 25 %, 10 %, 5 %, 2.5 % and 1.25 % with LEA SYMBOLS® and LEA NUMBERS®.

VA from 0.05 to 1.25 **(40 cm)**, Format: 15.8 cm x 15.8 cm

LEA SYMBOLS® – ring binder

With response chart and cover plate

Contrast: 2,5 %,

VA from 0.10 to 2.50 (3 m), Format: 12.7 cm x 12.7 cm

LEA SYMBOLS® – ring binder – crowding edging

Contrast: 2,5 %

VA from 0.20 to 1.00 (3 m), Format: $10.2 \text{ cm} \times 10.2 \text{ cm}$

LEA SYMBOLS® – Low Contrast ring binder Stock code: 52028

LEA NUMBERS® – Low Contrast ring binder

Random order, Format: 8.25 cm x 22.9 cm

SLOAN Low Contrast ring binder – letters

5 pages in 25 %, 10 %, 5 %, 2.5 % and 1.25 % contrast (1 \mathbf{m})

Format: 8.25 cm x 22.9 cm

LEA SYMBOLS® – Low Contrast ring binder Stock code: 52028A

(GL 251200)

LEA NUMBERS® – Low Contrast ring binder Stock code: 52075A

(GL 270500)

SLOAN Low Contrast ring binder – letters Stock code: 52182

5 pages in 25 %, 10 %, 5 %, 2.5 % and 1.25 % contrast **(1 m)** (GL 737000)

Stock code: 52222

Stock code: 52235

Stock code: 52236

Stock code: 52239

(GL 587700)

(GL 610800)

(GL 250625)

(GL 259225)

(GL 251100)

(GL 270400)

(GL 736000)

Stock code: 52075

Stock code: 52181



Accommodometer

Accommodometer German text Accommodometer English text

RAF accommodometer for assessment of accommodation and convergence near point.

Stock code: 53527 Stock code: 53527-E

Stock code: 53005

Stock code: 53019

Stock code: 53003

Stock code: 59520

Stock code: 53004



Bagolini hand-held frame and lenses

Bagolini hand-held frame (OCULUS®)

Cylinders in 45°/135° for testing of binocular vision under normal viewing conditions

Bagolini lenses (OCULUS®)

For trial frames (2 pieces), format: Ø 38 mm



Red glasses



For diplopia-provocation and angle-determination with the Maddox

cross.

Size: Ø 38 mm

Trapezoid dark Red glass

For tangent screen and examinations requiring a large glass diameter.

Size: 10 cm x 6 cm

Maddox cylinder (OCULUS®)

For use in trial frames, for measuring deviation at fusion break point.

Size: Ø 38 mm



#53001 #53002 #53016

Red-Green glasses, frames & flippers

Red-Green glasses (OCULUS®) Stock code: 53001

Red-Green hand-held frame (OCULUS®) Stock code: 53002

Red-Green flipper Stock code: 53016

4 lenses (GL 494300)

#53025

Wrap-around Red-Green glasses for adults Stock code: 53025

(GL 153160)

Wrap-around Red-Green glasses for children Stock code: 53026

(GL 153170)



Red-Green glasses for children with horse design Stock code: 53007 (GL 461300)

Red-Green glasses for children with butterfly design Stock code: 59538

Red-Green glasses for children with parrot design Stock code: 59539

(GL 461400)

(GL 461200)

Red-Green glasses for children with tiger design Stock code: 59540

(GL 461500)



Red-Green glasses with elastic strap

Stock code: 53014

Stock code: 53013

(GL 153090)




Red-Green tests



#53017

LED Red-Green test (Berens) Simple test to check binocular vision and fusion with three objects (dot/circle, figure, elephant).

Flashlight LEA SYMBOLS® Stock code: 53017



Stock code: 53018 **WORTH 4-Dot Flashlight**

(GL 466800)

(GL 520600)



Batteries required but not included For matching red-green glasses see above.

55

#59501 #59517 #59590 #59581 #59582

Occluders and Paddles

Classic Cover-Test Paddle (OCULUS®)

With clip

Spielmann Occluder Stock code: 59517

semi-transparent

Semi-transparent Occluder Stock code: 59590

high transparency (GL 760000)

Pinhole Occluder Stock code: 59581

1 central pinhole (Ø 1.5 mm) (GL 479500)

Pinhole Occluder Stock code: 59582

49 pinholes, long handle (GL 478500)

The images of the following items might differ slightly from the original.

Occluder red glass Stock code: 59507 (GL 206001)

Occluder Maddox cylinder Stock code: 59508

(GL 205001)

Stock code: 59501

Occluder Maddox cylinder / black Stock code: 59509

(GL 207000)

Occluder Maddox cylinder / red Stock code: 59510

(GL 209000)

Maddox Cylinder - central Pinhole Stock code: 59584

(GL 543500)

Occluders and Paddles



LEA™ Face disposable Occluders

50 pcs., can be given to the child following a health screening or a preventive paediatric check-up

Stock code: 59544

(GL 253200)



Plastic Occluder set, curved, grey

6 pcs.

(GL 750000)

Plastic Occluder set, curved, black

Stock code: 59550

Stock code: 59549

(GL 751000)



Plastic Occluder curved, grey

Stock code: 59547

(GL 700544)

Plastic Occluder curved, black

1 pc.

Stock code: 59548

(GL 700545)



Multiple Pinhole Spectacle Occluder

The holder enables rapid, complete coverage of the eye and visual acuity testing via multiple pinholes



Spectacle Occluder

The holder enables rapid, complete coverage of the eye

Stock code: 59546

Stock code: 59545

(GL 700012)

Polarisation/Red-green Flipper

4 lenses

Stock code: 53015

(GL 488400)



Clip-on Occluder, beige

To be mounted on the glasses frame (Set of 6 pieces)

Stock code: 53531

(GL 481150)



#53015

Clip-on Occluder, black

To be mounted on the glasses frame (1 piece)

Stock code: 59567

Stock code: 53536

(GL 121002)

(GL 166000)



Visual Field Fatigue Diffuser

For testing of VA and to conduct perimetry by fogging the contralateral eye. Plastic material (3 pieces)





Opaque Glass Occluders



Parrot opaque glass occluder

Stock code: 59536

(GL 461011)

Stock code: 59537

(GL 461111)

Tiger opaque glass occluder

Set of 2 frames (frosted lens) To enable visual acuity assessment even in the presence of nystagmus



Sunflower occlusion glasses

Both sides can flip up for VA testing.

Stock code: 59543

(GL 548500)



Flip-Up Occluder Glasses (frosted lens)

Flip-Up Occluder Glasses (opaque)

Both sides can flip up for VA testing. Size: 12.5 cm x 5 cm (width x height)

Lens diameter: 4.5 cm, temporals 11.8 cm length

Stock code: 59598

(GL 587611)

Stock code: 59599

(GL 587600)



#59542

Polarisation Occluder - butterfly

Polarisation Occluder - tiger

Polarisation Occluder - horse

Polarisation Occluder - parrot

For binocular and stereoscopic testing.

Stock code: 59533

(GL 460400)

Stock code: 59534

(GL 460700)

Stock code: 59541

(GL 462000)

Stock code: 59542

(GL 462100)



#52233

Paediatric Visual Testing Set



LEA SYMBOLS® / LEA NUMBERS® card #52018, LEA SYMBOLS®/ Tumbling E wallchart #52035, Butterfly opaque glass occluder #53009

Suitable stereo tests can be found on pages 68 and 69.



#52234


Provikit eye tests (Project Provision)

Details and composition on request (please refer to the illustration for the possible selection)

Stock code: 52234

Stock code: 52233

Visual Acuity Tests for patients with dementia





Please visit our online shop: www.eyesfirst.eu

OKN Instruments



The slow phase of the nystagmus is in the same direction as the movement of the flag or the drum. The nystagmus is named after the direction of the quick phase.

OKN flag in red-white

Size: 9.5 cm x 71 cm (GL 547400)

Stock code: 53502

Stock code: 53562

Stock code: 53510

Stock code: 53557

OKN flag in red-white

Size: 25 cm x 68 cm (GL 588600)

#53562



OKN drum for provocation and testing opto-kinetic nystagmus

(GL 547500) Height of drum: 25 cm

Overall height: 43 cm Diameter: 15 cm

OKN drum for provocation and testing opto-kinetic nystagmus

Height of drum: 10 cm Overall height: 28 cm

Diameter: 15 cm

(GL 616700)



Nystagmus glasses Frenzel

Nystagmus glasses (Frenzel) (DEHAG)

High plus lenses disrupt fixation and therefore enable detection of spontaneous nystagmus. The +15 D glasses provide magnification for easy assessment of eye movements. Battery version with illumination


Batteries required but not included



Maddox-Wing test

Maddox-Wing Test

Quick and efficient near test for horizontal, vertical and cyclodeviations. Differentiation between phoria and tropia must be carried out separately.



Maddox Cross

Maddox Cross with LED lights (OCULUS®)

Scale for 5 m and 1 m, incl. LED lights and transformer. For easy testing of horizontal and vertical deviations. Size: approx. 1.30 m (height / width)

Stock code: 59521

Stock code: 59519

Stock code: 53516

Trial Frames adults and children



#57019

Universal Trial frame – our "special brand"

Cost effective classic standard trial frame with a maximum of 5 lenses in either aperture. Ø 38 mm, PD 50 mm to 80 mm.

(GL 455500)

Stock code: 57019



Trial frame UTF-5080 PRO - our own "Premium" brand

Lightweight trial frame in a modern design, can accommodate 5 trial lenses (Ø 38 mm) per side (IPD 50 mm to 80 mm)

Stock code: 57057



#57001

Trial frame UB 3+ (OCULUS®)

Trial frame for clinics and medical practices, can hold a maximum of 5 lenses (Ø 38 mm) in either aperture (PD 48 mm to 80 mm).

Stock code: 57001



Trial frame UB4 (OCULUS®)

Modern trial frame in appealing design. R/L adjustable PD (48 to 80 mm). For up to 5 lenses (Ø 38 mm) in either aperture. Stock code: 57005



Trial frame UB 5 (OCULUS®)

Light-weight trial frame in attractive design. With removable coloured knobs and adjustable PD (40 to 72 mm).

For up to 5 small lenses (Ø 28 mm) in either aperture.

Stock code: 57006



#57006

Trial frame UB 6 (OCULUS®)

The new development focused on ergonomics, simple handling and optimal wear comfort. Can be fitted with 6 trial lenses (Ø 38 mm) in either aperture (PD 46 mm to 80 mm).

Stock code: 57049



Trial frame (OCULUS®)

Trial frame, robust

Light weight trial frame for up to 3 lenses (Ø 38 mm) in either aperture. For refraction with adjustable PD (59 mm to 67 mm) and nose-pad.

Stock code: 57053



#57046

Trial frame in robust quality (Ø 38 mm) with adjustable PD (54 mm to 70 mm).

Stock code: 57046

(GL 469900)

#57046-A



Trial Frames children

Baby trial frame COMO (PD 50 mm)

For up to 2 lenses (Ø 38 mm) in either aperture.

Toddler trial frame COMO (PD 57 mm)

For up to 2 lenses (Ø 38 mm) in either aperture.

Children trial frame, robust (image similar)

Trial frame for children in robust quality (Ø 38 mm) with adjustable PD 48 mm to 60 mm.

Stock code: 57050

Stock code: 57052 (GL 472400)

Stock code: 57046-A

Stock code: 57004

Stock code: 57015A

Stock code: 57015B

(GL 472300)

Children trial frame (OCULUS®)

For up to 3 lenses (Ø 38 mm) in either aperture. For refraction with adjustable PD 54 mm to 58 mm and nose-pad.



Refraction Trial Frame sets

Refraction set for children

4 trial frames (Ø 38 mm) **with fixed PD** in the following sizes: PD 50, 52, 54, 56 mm

Refraction set for adolescents and adults

6 trial frames (Ø 38 mm) **with fixed PD** in the following sizes: PD 58, 60, 62, 64, 66, 68 mm

Refraction Glasses children









Useful for retinoscopy, testing of convergence excess and AC:A ratio.

Refraction glasses tiger	+1.50 D	(GL 528900)	Stock code: 57010
Refraction glasses butterfly	+1.50 D	(GL 516000)	Stock code: 57039
Refraction glasses horse	+1.50 D	(GL 525900)	Stock code: 57041
Refraction glasses parrot	+1.50 D	(GL 525700)	Stock code: 57040
Refraction glasses tiger	+2.00 D	(GL 491300)	Stock code: 57011
Refraction glasses butterfly	+2.00 D	(GL 491000)	Stock code: 57029
Refraction glasses horse	+2.00 D	(GL 491100)	Stock code: 57030
Refraction glasses parrot	+2.00 D	(GL 491200)	Stock code: 57031
Refraction glasses tiger	+2.50 D	(GL 500611)	Stock code: 57035
Refraction glasses butterfly	+2.50 D	(GL 500300)	Stock code: 57032
Refraction glasses horse	+2.50 D	(GL 500400)	Stock code: 57033
Refraction glasses parrot	+2.50 D	(GL 500500)	Stock code: 57034
Refraction glasses tiger	+3.00 D	(GL 515200)	Stock code: 57012
Refraction glasses butterfly	+3.00 D	(GL 514900)	Stock code: 57036
Refraction glasses horse	+3.00 D	(GL 515100)	Stock code: 57037
Refraction glasses parrot	+3.00 D	(GL 515000)	Stock code: 57038



Trial Lenses set

Trial Lenses set

90 lenses (Ø 38 mm) with metal rim, spherical 0.25 D to 10.00 D, minus cylinders 0.25 D to 3.00 D, 6 additional lenses (Occluder, Pinhole, Slit, Red Filter, Maddox Red Filter, Green Filter) plus Jackson cross cylinder 0.25 D with handle (all lenses without anti reflective coating).

For a selection of single lenses and additional lenses see page 66

Jackson Cross cylinders



For subjective refraction

Jackson cross cylinder ± 0.25 D Jackson cross cylinder ± 0.50 D Jackson cross cylinder ± 0.75 D Jackson cross cylinder ± 1.00 D

ross cylinder ±0.75 D Stock code: 57043 ross cylinder ±1.00 D Stock code: 57044

Confirmation bar



Confirmation bar (OCULUS®)


For spherical fine adjustment in subjective refraction. ± 0.25 D und ± 0.50 D

Retinoscopy Bar sets



Retinoscopy bar set (Luneau)

Retinoscopy bars in parallel arrangement, one plus and minus bar each from 0.5 D to 15.0 D, incl. case, Powers per bar: 0.5; 1.0; 1.5; 2.0; 2.5; 3.0; 3.5; 4.0; 5.0; 6.0; 7.0; 8.0; 9.0; 10.0; 12.0; 15.0 D Glass diameter approx. 15 mm



Powers per bar: 0.5; 1.0; 1.5; 2.0; 2.5; 3.0; 3.5; 4.0; 4.5; 5.0; 5.5; 6.0; 6.5; 7.0; 8.0 and 12 D.

One plus and minus bar of 0.5 D to 12.0 D each, in case. The 180° Tabo notation above the lenses facilitates determination of the cylinder axis position.

Glass diameter approx. 16 mm

Stock code: 57020

Stock code: 57000

Stock code: 57045

Stock code: 57002

Stock code: 57003

Stock code: 57008









Confirmation Flippers spherical

Binocular confirmation flippers of various power, each with 2 parallel plus and minus lenses

Confirmation flipper ± 0.25 D	(GL 127025)	Stock code: 57021
Confirmation flipper ± 0.50 D	(GL 127050)	Stock code: 57022
Confirmation flipper ± 0.75 D	(GL 127075)	Stock code: 57023
Confirmation flipper ±1.00 D	(GL 127100)	Stock code: 57024
Confirmation flipper ±1.25 D	(GL 127125)	Stock code: 57025
Confirmation flipper ±1.50 D	(GL 127150)	Stock code: 57026
Confirmation flipper ±1.75 D	(GL 127175)	Stock code: 57027
Confirmation flipper ± 2.00 D	(GL 127200)	Stock code: 57028
Confirmation flipper ±2.50 D	(GL 127250)	Stock code: 57047
Confirmation flipper ± 3.00 D	(GL 127300)	Stock code: 57048

Confirmation flipper set including stand

Stock code: 57051

Stock code: 57054

(GL 499700)

(GL 552500)

Confirmation flipper set in carrying case (not pictured)

6 pieces, powers: 0.25 D, 0.50 D, 0.75 D, 1.00 D, 1.50 D, 2.00 D



Prism Flippers

Prism flipper 1△	(GL 137001)	Stock code: 55435
Prism flipper 2△	(GL 137002)	Stock code: 55429
Prism flipper 3△	(GL 137003)	Stock code: 55430
Prism flipper 4△	(GL 137004)	Stock code: 55431
Prism flipper 6△	(GL 137006)	Stock code: 55432
Prism flipper 8△	(GL 137008)	Stock code: 55433

The indicated power relates to the single prism.



Prism Flipper Set (image similar) Set with 6 prism powers:

Powers: $1\triangle$, $2\triangle$, $3\triangle$, $4\triangle$, $6\triangle$, $8\triangle$

Stand for 6 Confirmation/Prism Flippers Stock code: 59964

(GL 501600)

(GL 624500)

Stock code: 55434

REFRACTION $GL = GOOD\text{-}LITE^{\otimes}$ stock code



Metal Rim lenses

Suitable for perimetry and any examinations where a large visual field is required (38 mm diameter).

	Stock code:		Stock code:		Stock code:
Metal Rim lens	+0.25 D 57200	Metal Rim lens, cylinder	−3.75 D 57244	Metal Rim lens	+13.00 D 57288
Metal Rim lens	-0.25 D 57201	Metal Rim lens	+4.00 D 57245	Metal Rim lens	-13.00 D 57289
Metal Rim lens, cylinder	-0.25 D 57202	Metal Rim lens	-4.00 D 57246	Metal Rim lens	+14.00 D 57290
Metal Rim lens	+0.50 D 57203	Metal Rim lens, cylinder	-4.00 D 57247	Metal Rim lens	-14.00 D 57291
Metal Rim lens	−0.50 D 57204	Metal Rim lens	+4.25 D 57248	Metal Rim lens	+15.00 D 57292
Metal Rim lens, cylinder	−0.50 D 57205	Metal Rim lens	-4.25 D 57249	Metal Rim lens	-15.00 D 57293
Metal Rim lens	+0.75 D 57206	Metal Rim lens	+4.50 D 57250	Metal Rim lens	+16.00 D 57294
Metal Rim lens	-0.75 D 57207	Metal Rim lens	-4.50 D 57251	Metal Rim lens	-16.00 D 57295
Metal Rim lens, cylinder	-0.75 D 57208	Metal Rim lens, Cylinder	-4.50 D 57252	Metal Rim lens	+18.00 D 57296
Metal Rim lens	+1.00 D 57209	Metal Rim lens	+4.75 D 57253	Metal Rim lens	–18.00 D 57297
Metal Rim lens	–1.00 D 57210	Metal Rim lens	-4.75 D 57254	Metal Rim lens	+20.00 D 57298
Metal Rim lens, cylinder	–1.00 D 57211	Metal Rim lens	+5.00 D 57255	Metal Rim lens	-20.00 D 57299
Metal Rim lens	+1.25 D 57212	Metal Rim lens	−5.00 D 57256		
Metal Rim lens	–1.25 D 57213	Metal Rim lens, cylinder	−5.00 D 57257	Plano lens	57300
Metal Rim lens, cylinder	–1.25 D 57214	Metal Rim lens	+5.25 D 57258	Frosted lens	57301
Metal Rim lens	+1.50 D 57215	Metal Rim lens	−5.50 D 57259	Occluder, black	57302
Metal Rim lens	–1.50 D 57216	Metal Rim lens	+5.50 D 57260	Red filter	57303
Metal Rim lens, cylinder	–1.50 D 57217	Metal Rim lens	−5.50 D 57261	Green filter	57304
Metal Rim lens	+1.75 D 57218	Metal Rim lens, cylinder	−5.50 D 57262	Pinhole	1.5 mm 57305
Metal Rim lens	–1.75 D 57219	Metal Rim lens	+5.75 D 57263	Pinhole	1.0 mm 57306
Metal Rim lens, cylinder	–1.75 D 57220	Metal Rim lens	-5.75 D 57264	Slit	1.5 mm 57307
Metal Rim lens	+2.00 D 57221	Metal Rim lens	+6.00 D 57265	Slit	1.0 mm 57308
Metal Rim lens	-2.00 D 57222	Metal Rim lens	−6.00 D 57266	Marking cross	57310
Metal Rim lens, cylinder	−2.00 D 57223	Metal Rim lens, cylinder	−6.00 D 57267		
Metal Rim lens	+2.25 D 57224	Metal Rim lens	+6.50 D 57268	Prism lens	0.50△ 57311
Metal Rim lens	–2.25 D 57225	Metal Rim lens	−6.50 D 57269	Prism lens	1.00△ 57312
Metal Rim lens, cylinder	–2.25 D 57226	Metal Rim lens	+7.00 D 57270	Prism lens	2.00 △ 57313
Metal Rim lens	+2.50 D 57227	Metal Rim lens	–7.00 D 57271	Prism lens	3.00△ 57314
Metal Rim lens	–2.50 D 57228	Metal Rim lens	+7.50 D 57272	Prism lens	4.00△ 57315
Metal Rim lens, cylinder	−2.50 D 57229	Metal Rim lens	–7.50 D 57273	Prism lens	5.00△ 57316
Metal Rim lens	+2.75 D 57230	Metal Rim lens	+8.00 D 57274	Prism lens	6.00△ 57317
Metal Rim lens	–2.75 D 57231	Metal Rim lens	−8.00 D 57275	Prism lens	8.00△ 57318
Metal Rim lens, cylinder	-2.75 D 57232	Metal Rim lens	+8.50 D 57276	Prism lens	10.00△ 57319
Metal Rim lens	+3.00 D 57233	Metal Rim lens	−8.50 D 57277		
Metal Rim lens	−3.00 D 57234	Metal Rim lens	+9.00 D 57278	Maddox Cylinder	white 57309
Metal Rim lens, cylinder	−3.00 D 57235	Metal Rim lens	−9.00 D 57279		(GL 473200)
Metal Rim lens	+3.25 D 57236	Metal Rim lens	+9.50 D 57280	Maddox Cylinder	red 57320
Metal Rim lens	−3.25 D 57237	Metal Rim lens	−9.50 D 57281		(GL 473000)
Metal Rim lens, cylinder	−3.25 D 57238	Metal Rim lens	+10.00 D 57282		
Metal Rim lens	+3.50 D 57239	Metal Rim lens	-10.00 D 57283		
Metal Rim lens	−3.50 D 57240	Metal Rim lens	+11.00 D 57284		
Metal Rim lens, cylinder	−3.50 D 57241	Metal Rim lens	-11.00 D 57285		
Metal Rim lens	+3.75 D 57242	Metal Rim lens	+12.00 D 57286	All lenses have a di	amotor of 20 mm
Metal Rim lens	−3.75 D 57243	Metal Rim lens	-12.00 D 57287	All letises flave a al	umerer or so IIIIII.

Retina Lenses MediWorks



Retina lens 20 D

High degree of light transmission, double aspherical precision design, large observation field (material: glass).

Magnification: 3.13 x, Lens diameter: 48 mm, Visual field: 46°;

Working distance: 50 mm



Retina lens 78 D

High degree of light transmission, double aspherical precision design, large observation field (material: glass).

Magnification: 0.93 x, Lens diameter: 26 mm, Visual field: 81.4°

Working distance: 8 mm



Retina lens 90 D

High degree of light transmission, double aspherical precision design, large observation field (material: glass).

Magnification: 0.76 x, Lens diameter: 21 mm, Visual field: 74.2°

Working distance: 7 mm



#57702

Stock code: 57700

Stock code: 57701

Stock code: 57702



LANG Stereo tests

NEW!

LANG Test I-R new edition

Stock code: 54000-R

#54000-R

3 objects with 1200, 600 and 400, suitable for use with stereopsis screening from approx.

age 12 months.



NEW!

LANG Stereotest® II

Stock code: 54001 Stock code: 54001-R

LANG-Stereotest® II-R

4 objects with 600, 400 and 200, with a control object that is also detected monocularly. Suitable for use with stereopsis screening from

approx. age 12 months old.





#54021

LANG STEREOPAD®

Stock code: 54021

6 magnetic test cards from 1000" to 50", which can be arranged individually or in groups. Suitable for stereo screening from 12 months of age on a preferential looking basis.

Monopod - support for LANG-STEREOPAD®

Stock code: 54021-S

Stock code: 54003

Stock code: 54005



No glasses are needed to perform the LANG stereo tests.

Stereo tests diverse



Titmus stereo test fly (Stereo Optical) Fly, rings, animals, from 3552" to 40" Polarisation glasses included



Titmus stereo test butterfly (Stereo Optical) Butterfly, rings, animals from 2500" to 40" Polarisation glasses included



Randot stereo test (Stereo Optical) Symbols, rings, animals, from 500" to 20" Polarisation glasses included

Stock code: 54004





NEW!



Stereo tests diverse

Stereo test "Fly" with LEA SYMBOLS® (Stereo Optical)

Random dot Stereo test with fly, rings and LEA SYMBOLS®, from 4800" to 20" With legend on the back and polarisation glasses

Stereo test "Butterfly" with LEA SYMBOLS® (Stereo Optical)

Random dot Stereo test with butterfly, rings and LEA SYMBOLS $^{\! \otimes}\! ,$

from 2500" to 20"

With legend on the back and polarisation glasses

Dinosaur Randot stereo test (Bernell)

New, glasses free stereo test with lenticular technology. Images can either appear to protrude outwards or recede inwards.

Stereopsis from 400" to 40"

TNO Test Stock code: 54002

With Red-Green glasses, different tests for binocularity and stereopsis in a ring binder. Contains 6 plates

TNO Test without Red-Green glasses) Stock code: 54002-A

Availability upon request

Availability

upon request

Stock code: 54006

Stock code: 54007

Stock code: 54020

(GL 370129)

TNO Red-Green glasses (red = left) Stock code: 54002-B

Availability upon request

Please note: the TNO glasses have the red filter on the left side.

Colour tests

Ishihara colour test with 38 plates

The classic colour test for quick reference, with numbers and shapes.

Ishihara colour test with 24 plates

Ishihara colour test with 14 plates

Stock code: 56000

Stock code: 56013

Stock code: 56010

Stock code: 56001

(GL 730022)

Farnsworth D-15 (Luneau)

Abbreviated Farnsworth 100 Hue test as short screening test for differentiation between slight and distinct colour blindness.

Farnsworth D-15 (GOOD-LITE®) Stock code: 56019

Serves as a quick screening test for the differentiation of mild and advanced colour deficencies.

Lanthony colour test (GOOD-LITE®) Stock code: 56021

Desaturated Farnsworth D-15 test (GL 730023)

Velhagen plates Stock code: 56012

Book with 36 testing plates for quick testing of colour vision disorders. Also suitable for use with children and illiterate individuals. Also includes tests for blue-yellow colour vision disorders.



#56012 Tafeln zur Pr Farbensinne

Control of the contro



Colour tests

LEA™ Panel 16 – Colour test

With large colour discs (\emptyset 3.3 cm) for testing of Low-Vision patients Travel Tube with black cloth for testing. Includes rings to reduce aperture for general testing

Format Travel Tube: Ø 6 cm, Length: 22 cm




(GL 260260)



LEA™ Panel 16 - Colour test

With large colour discs (Ø 3.3 cm) for testing of Low-Vision patients, in a transparent box

Format Box: 74 cm x 8.5 cm





LEA™ Panel 16 - Colour test, double set

The double set makes it easier to perform the colour test, especially on children. In this test, the colour patterns are laid out in any sequence and the child is asked to allocate the matching colour values from the second set.

Format Box: 74 cm x 8.5 cm, colour disc diameter of 3.3 cm

Stock code: 56016 (GL 260100)



LEA™ Panel 16 – Colour test – double set in a Travel Tube

With large colour discs (Ø 3.3 cm) for testing of Low-Vision patients. Format: 2 Tubes à Ø 6 cm, Length: 22 cm

(GL 260160)



Waggoner colour test "Color Vision Screening Made Easy"


The popular colour test with symbols is now available again in a new version. Contains 1 demonstration card and 9 test cards with circle, star, square for rapid screening of colour vision and 4 additional charts for screening.

Stock code: 56015

Stock code: 56024

(GL 730005)

Stock code: 56017



HRR Colour test, laminated version (GOOD-LITE®)

The HRR (Hardy Rand and Rittler) Standard Pseudoisochromatic Test, 4th Edition provides several important features to provide the most advanced color vision test available: congenital and acquired testing, identification of the type of defect and diagnosis of the extent of the defect as well as quick positive classification of normals. The figures (circle, cross, triangle) used by the HRR Pseudoisochromatic Plates are independent of language and suitable for both adults and children. Contents: 4 demonstration plates, 6 plates for Screening, 14 plates to diagnose color deficiency.



HRR Colour Test / Amsler Test, laminated (GOOD-LITE®)

Ring binder with HRR Colour Test pages and Amsler Test pages.

Stock code: 56028

(GL 730020)



71

Visual Field tests

LEA FLICKER WAND®

Improves Confrontation Visual Field testing. With flickering or nonflickering stimulus. The curved wand is moved from behind the person forward in all four quadrants of VF to detect losses in the periphery of the VF. In particular useful in assessment of Low Vision and in Rehabilitation.

Format: 67.3 cm length in total.

Stock code: 53511

(GL 280000)

Damato Campimeter

Hand-held, portable visual field test, reciprocally printed. Ideal for testing in cases where normal VF testing is impractical. Tests 30° central VF by presenting a black-white stimulus. The patient is asked to fixate the relevant dots marked by numbers.

Format: 41 cm x 41 cm

Stock code: 53508

(GL 293000)

#53508



#59812

Amsler test ring binder (GOOD-LITE®)

7 pages: Standard (Amsler Chart), Central (Scotoma Diagonal), Color Scotoma (Red on Black), No Lines (Dots), Metamorphopsia (Parallel Lines on Black), Metamorphopsia Reading card (Added Central Lines), Juxta-central (Standard with Small Squares).

Stock code: 59812

(GL 730026)



HRR Colour Test / Amsler Test, laminated (GOOD-LITE®) Ring binder with HRR Colour Test pages and Amsler Test pages. Stock code: 56028 (GL 730020)

Consumables



Single-use Tonometer Prisms L900 (Luneau)

For all types of Goldmann applanation tonometers incl. Perkins tonometer. The individually sterile-packaged tonometer heads prevent cross contamination and are particularly useful in infectious diseases and after surgeries. Excellent image quality, simple handling. 100 pieces/box

Stock code: 53700

Tiered prices for large quantities available on request

Application graphic for download





Visual Field Fatigue Diffuser

For testing of VA and to conduct perimetry by fogging the contralateral eye. Plastic material (3 pieces)

Stock code: 53536

(GL 166000)



Pro-ophta® eye swabs

500 pieces (64 mm x 5 mm), non-sterile

Stock code: 58008



#53540

Aesthesiometer Cochet-Bonnet (Luneau)

Classic Aesthesiometer to measure corneal sensitivity with nylon thread 12/100 mm.

Replacement fil for Aesthesiometer (Luneau)

Caliper: 12/100 mm

Stock code: 53540-E

Stock code: 59596

Stock code: 59602

(GL 108100)

Stock code: 53540



Marsden Fixation Ball

3.2" hollow ball made of durable vinyl with applied black letters and numbers on a 10 foot cord.

For testing and training of eye movements.



QuikVue® Plus VPA-200

smartphone adapter including yellow filter for anterior segment photography with 10x magnification. Three LED illumination modes to choose from: White in two brightness levels, blue light including switchable yellow filter. Li-ion battery with charging option via micro USB cable. Ideal for mobile image documentation, telemedicine and patient education. Weight: approx. 59 g



Before placing your order, please confirm whether your smartphone has the appropriate camera position for the QuikVue® adapter.

Further information about the download





Multiple Pinhole Spectacle Occluder

Stock code: 59545

(GL 700010)

Stock code: 59546

(GL 700012)

Spectacle Occluder without pinhol	Spectac	le Occ	luder	with	out (pinho	le
-----------------------------------	---------	--------	-------	------	-------	-------	----

The holder enables rapid, complete coverage of an eye for visual acuity testing via multiple pinholes (without image).

#59506 #59528 #59514

Polarisation Glasses for adults

suitable in combination with Stereo tests

Polarisation Glasses, intermediate size Stock code: 59528

suitable in combination with Stereo tests

Polarisation Glasses for children

suitable in combination with Stereo tests

Stock code: 59514

Stock code: 57016

Stock code: 57018

Stock code: 57059

Stock code: 59506

(GL 103050)

(GL 103175)

(GL 103150)

#57016



Nose pad blue, soft (OCULUS®)

3 pieces, for OCULUS® trial frames

Nose pad black, soft (OCULUS®)

2 pieces, for classic OCULUS® UB 3 only

Nose pad set grey grau (OCULUS®)

Stock code: 57017 (OC 42183)

(OC 42511)

#57018



Nose pad for OCULUS® children's trial frame (OCULUS®)

1 piece

(OC 42483)

set of 3, for OCULUS® trial frame UB 4 (OC 42510)

75



NEW! #53515

LEA™ book "What and How Does This Child See? (2nd edition) The book covers different aspects of visual function and maturation (in English).

Stock code: 53515

(GL 791000)

LEA™ USB stick "What and How Does This Child See?

USB version of the book (in English)

Stock code: 53542

Stock code: 53517

(GL 792000)

(GL 790000)

#53517



#52032

#52033

#52242

#52034

#52031

101 101

101 101

#52241

DVD "Leo Learns By Doing"

LEA™ Symbols response plate

4 objects/plate, semi-transparent

LEA™ Symbols response cards

LEA™ Numbers response plate

LEA™ Symbols response plate – crowding edging

4 single cards

4 objects/plate

4 Numbers/plate

4 single cards

4 objects/plate

4 objects/plate

HOTV response cards

HOTV response plate

The DVD tells the story of a visually impaired child named Leo. Step by step, the therapist explains techniques that parents can use to develop and support Leo's existing abilities. Soundtrack in English, Finnish, Swedish, Russian, Chinese and Spanish.

Stock code: 52031

(GL 251800)

(GL 255303)

(GL 270030)

(GL 700515)

Stock code: 52034

Stock code: 52245

(GL 251701)

Stock code: 52032

Stock code: 52242

Stock code: 52241

Stock code: 52033

(GL 700523)

(GL 700530)

Stock code: 57013



Reversible Occluding Glasses

HOTV response plate - crowding edging

Occluding glasses for rapid reversible occlusion in vision tests and monocular exercises.

(GL 452300)



Hertel Mirror Exophthalmometer (OCULUS®)

Classic mirror exophthalmometer for measuring protrusion in bilateral comparison. Scale with mm increments

Stock code: 53525



Neutral Density filter (Gulden)

To examine aff erent pupillary defects (APD). Also suitable as alternative to the Sbisa bar.

Filter values: 0.3; 0.6; 0.9; 1.2; 1.5; 1.8 density log units.

Stock code: 53010 (GL 151170)

Pupillometer (Bernell)

The front shows different pupil diameters, the back lists expected pupil size in different lighting.

Set of 5 pieces

Format: Ø approx. 83 mm

Stock code: 53521 (GL 566500)

#53535

Flashing Fixation Star for Retinoscopes/Ophthamoscopes

Set of 3 for Welch Alynn® and Keeler Ophthalmoscopes and Retinoscopes. Comes with double-sided tape for mounting.

Stock code: 53535

(GL 489000)

#59588

#53521



Diagnostic Torch Penlight

Handy diagnostic torch with replaceable batteries.

Stock code: 59588

#59959



Brock string with 5 beads

Package of 12 strings, length 3 m.

Single Brock string

Stock code: 59959

(GL 593800)

Stock code: 59959-A

(GL 603200)



Indicator - E-hook (OCULUS®)

Stock code: 53500

Indicator - Landolt ring (OCULUS®)

Stock code: 53501

Landolt-C Response Optotype

Made of durable plastic Format: Ø approx. 15 cm Stock code: 59577

(GL 602600)

#59585



Tumbling E Response Optotype

Made of durable plastic Format: approx. 15 cm x 15 cm Stock code: 59585

(GL 673700)

LEA™ Ruler Pads

Pad of 25 sheets as give away for patients. Helps to practice the LEA symbols response at home.

Stock code: 149

Stock code: 59968

QR code for multi-lingual product information regarding LEA™ tests





Eye model, 6-part

Eye model standing on base, 3 times full-size. Separates into two halves with retina, iris, lens and vitreous body.

nerve in natural position in the eye socket. Separates into two halves

with cornea, muscle insertions, retina, iris, lens and vitreous body.

Eye model in eye socket, 7-part Stock code: 59969 Eye model in eye socket, 3 times full-size. Presentation with optic



78

The ORTOPAD selection







Ophthalmoscopes



● • . | •

HEINE BETA® 200S Ophthalmoscope*

74 different diopter-steps, aspherical optics. 7 different filter, incl. fixation-star with polar coordinate, 3.5 V.

Stock code: 57500

(HEINE C-002.30.120)



HEINE BETA® 200S LED Ophthalmoscope*

74 different diopter-steps, aspherical optics. 7 different filter, incl. fixation-star with polar coordinate, 3.5 V. LED Ophthalmoscope with continuous brightness control between 100 % and 3 %.

Stock code: 57500-L (HEINE C-008.30.120)



HEINE BETA® 200 Ophthalmoscope*

28 different diopter-steps, aspherical optics, 6 different filter, incl. fixation star with polar coordinate, 3.5 V. (image shows ophthalmoscope with separate handle)

Stock code: 57501 (HEINE C-002.30.100)



HEINE BETA® 200 LED Ophthalmoscope*

28 different diopter-steps, aspherical optics, 6 different filter, incl. fixation star with polar coordinate, 3.5 V. LED Ophthalmoscope with continuous brightness control between 100 % and 3 %.

Stock code: 57501-L (HEINE C-008.30.100)



HEINE BETA® 200 Ophthalmoscope M2 filters*

27 different diopter-steps, aspherical optics, 6 different filter, incl. fixation star without polar coordinate.

Stock code: 57501A

(HEINE C-002.30.102)



HEINE K 180® Ophthalmoskop*

27 different diopter-steps, spherical optics, as standard-Ophthalmoscope. 5 different filter, incl. fixation star with polar coordinate, 3.5 V.

Stock code: 57502

(HEINE C-002.30.206)

#57556



HEINE mini 3000® LED Ophthalmoscope

Handy ophthalmoscope with LED lighting, and batteries housed in handle (standard AA size). Maintenance-free, since no lamp replacement is necessary.

Stock code: 57556 (HEINE D-008.71.120)

#57541



HEINE® Indirect Ophthalmoscope, handheld*

Monocular ophthalmoscope in pocket format for the ward round and when out and about. One-handed operation with integrated interference red-free filter and + 3 D loupe attachment. Additional binocular attachment available. Fitting 3.5 V rechargeable handles incl. Lens attachment.

Stock code: 57541 (HEINE C-002.33.001)

For the matching handles and accessories please refer to page 83.



Ophthalmology Sets

#57570



Diagnostic Set HEINE BETA® 200 LED

BETA® 200 Ophthalmoscope with LED, Retinoscope (LED), BETA 4 USB rechargeable handle and USB cord with plug-in power supply. Including hard box

Stock code: 57570

(HEINE C-145.28.388)

Stock code: 57567

(HEINE C-182.10.118)

Stock code: 57503-L

HEINE K 180® Ophthalmology set

Ophthalmoscope with battery handle, 2.5 $\rm V$

This handle does not fit the chargers and uses standard batteries size "C".

(batteries not included)

#57567



Retinoscopes

HEINE BETA® 200 LED Retinoscope*

LED retinoscope with a brilliant streak image. Completely homogeneous (HEINE C-008.15.353) and bright light for a precise, fast diagnosis.



For the matching handles and accessories please refer to page 84 to 85.

#57503-L

* without handle HEINE 8



Hand-held Instruments

HEINE® Lambda 100 Retinometer*

For determining potential visual acuity in the case of cataracts and other opacities. Can be used anywhere, even with bedridden patients. First compact instrument of this type. XHL xenon halogen technology with brightness control.

With decimal scale 3.5 V

Stock code: 57563 (HEINE C-002.35.010)

With fraction scale 3.5 V

Stock code: 57564

Patient card HEINE® Lambda 100 Retinometer

for optimal compliance

#57565

(HEINE C-002.35.015) Stock code: 57565

(HEINE C-000.35.005)

Stock code: 57552

(HEINE C-002.14.602)

HEINE HSL 150® hand-held Slit Lamp

Monocular hand-held Slit Lamp for examination of the anterior eye segment. Slit size infinitely adjustable, 6 × magnification, blue interference filter.

With Xenon-halogen lamp. Fits all BETA® handles with 3.5 V. The headpiece HSL 150 weighs only approx. 70 g; the lightest slit lamp of its type.



#57552

#57543

HEINE® Ophthalmic Examination Lamp*

For examination of the anterior segment in focal light. Five different luminous-field diaphragms, suitable for detecting erosions, oedemas and foreign bodies, 3.5 V. Fits all BETA® handles.

Stock code: 57543

(HEINE C-002.14.400)



HEINE mini-c® Clip Lamp LED

Long-life pocket diagnostic light in compact, modern design with chrome head. Bright, focused light.

Stock code: 57539

(HEINE D-008.73.109)

HEINE® Finoff Transilluminator*

For scleral transillumination, 3.5 V.

Stock code: 57568 (HEINE C-002.17.080)



HEINE® Charger NT 4

Charger for BETA® 4 NT handles (not suitable for old type NT rechargeable handles and USB rechargeable handle). (Heine X-002.99.494)

HEINE® Charger NT 300®

Charger incl. adapter for BETA® handles. LED rings indicate charging process.

Chargers



Stock code: 57537 (HEINE X-002.99.495)

Stock code: 57569

(HEINE X-002.99.494)

HEINE



Handles

Handles for chargers:



#57573



#5750*6*

Pictures #57573 and #57506 show the underside of the respective charging handle.



HEINE® Rechargeable handle BETA® 4 NT

For NT 4 charger, including Li-lon rechargeable battery and base unit BETA® 4 NT.

HEINE® Rechargeable handle BETA® NT

Suitable for the NT 300® charger (not for NT 4) with NiMH battery 3.5 V

HEINE® Rechargeable handle BETA® L

Suitable for the NT 300® charger (not for NT 4) with Li-Ion battery 3.5 V

Stock code: 57573 (HEINE X-007.99.396)

Stock code: 57505 (HEINE X-002.99.411)

Stock code: 57506

510CK Code: 5/506

(HEINE X-007.99.395)

Specifically suitable for portable use:



Rechargeable handle including Li-Ion rechargeable battery, BETA® 4 USB base unit, USB cord, plug-in power supply and handle rest. Does not need a charger.

HEINE® Battery handle

For 2 alkali-batteries, size "C", 2.5 V; **not** suitable for BETA® LED ophthalmoscopes

HEINE® Large battery handle

For 2 alkali-batteries, size "D", 2.5 V; **not** suitable for BETA® LED ophthalmoscopes (without batteries)

Stock code: 57571

(HEINE X-007.99.388)

Stock code: 57531 (HEINE X-001.99.118)

Stock code: 57509

(HEINE X-001.99.120)

Existing BETA® handles can be upgraded with an additional base unit to the BETA® 4 USB or BETA® 4 NT charging handles.

#57509

Rechargeable Batteries

#57538



#57512



HEINE® Li-Ion/Li-Ion L rechargeable battery
For RETA® L handles with integrated charging

For BETA® L handles with integrated charging status display in the base unit and rechargeable handle BETA 4 NT and BETA 4 USB (**not** suitable for charging handles BETA® NT, BETA® R and BETA® TR)

HEINE® NiMH rechargeable battery 3.5 V

For BETA® NT, BETA® R and BETA® TR handles

Forms and and bindly an acify the time of the battery (a.e. V. 00 00

Stock code: 57538 (HEINE X-007.99.383)

Stock code: 57512 (HEINE X-002.99.315)

For re-orders kindly specify the type of the battery (e.g.: X-02.99...).



Spare bulbs

		For HEINE BETA® 200S, BETA® 200 and BETA® M2	
		Ophthalmoscope	
		2.5 V	Stock code: 57518 (HEINE X-001.88.069)
#57519	HEINE XHLV	3.5 V	Stock code: 57519 (HEINE X-002.88.070)
		For HEINE K180® Ophthalmoscope	
		2.5 V	Stock code: 57520 (HEINE X-001.88.084)
#57521	HEINE XHI #086 3,5V	3.5 V	Stock code: 57521 (HEINE X-002.88.086)
		For HEINE BETA 200® streak Retinoscope	
		2.5 V	Stock code: 57522 (HEINE X-001.88.087)
#57523	(055 53)	3.5 V	Stock code: 57523 (HEINE X-002.88.089)
		For HEINE BETA 200® spot Retinoscope	
		2.5 V	Stock code: 57524 (HEINE X-001.88.088)
#57528		For HEINE® Diagnostic Torch ClipLight® 2.5 V	Stock code: 57528 (HEINE X-001.88.094)
#57540	(later Cas)	For HEINE® mini-c® Clip Lamp Xenon-Halogen-Lampe, 2.5 V	Stock code: 57540 (HEINE X-001.88.108)
#57548	HEINE XHI #050 3,5V	For HEINE® Indirect Ophthalmoscope 3.5 V	Stock code: 57548 (HEINE X-002.88.050)
#57553	(MONE SAIL)	For HEINE® hand-held Slit Lamp HSL 150® 3.5 V	Stock code: 57553 (HEINE X-002.88.099)
57529	HEINE XHL 1048 9,59	For HEINE® NIKATRON Retinoscope 3.5 V	Stock code: 57529 (HEINE X-002.88.048)
#57536	9647 3.5V	For HEINE® Ophthalmic Examination Lamp and HEINE® AUTOFOC Ophthalmoscope $3.5\ \lor$	Stock code: 57536 (HEINE X-002.88.047)
#57566	HEINE XHL #078 3,5V	For HEINE® Lambda Retinometer and HEINE® Finoff Transilluminator	Stock code: 57566 (HEINE X-002.88.078)

For re-orders kindly specify the respective no. of the bulb as indicated on the bulb (e.g.: #089)

XHL Xenon Halogen (XHL®), 3.5 V





Accessories

Loupe attachment 10 x for hand-held Slit Lamp HSL 150®

Stock code: 57555

(HEINE C-000.14.606)

Blue filter for hand-held Slit Lamp HSL 150®

Stock code: 57554 (HEINE C-000.14.605)

Stock code: 57546

Stock code: 57532

Stock code: 57535

Stock code: 57565

Stock code: 57547

Stock code: 57544

Stock code: 57545 (HEINE C-000.33.011)

Fundus charts

(HEINE C-000.33.208)

Block of 50 pages

HEINE® Fixation Cards with holder

(HEINE C-000.15.360)

For dynamic retinoscopy to be fixed on the BETA® 200 retinoscope.

Orange filter

(HEINE C-000.15.359) For BETA 200® retinoscope to reduce reflections

Patient card HEINE® Lambda 100 Retinometer

(HEINE C-000.35.005) for optimal compliance

Hardbox for Ophthalmology sets

BETA® handles and BETA® R handles fit the box

For ophtalmoscope, retinoscope, rechargeable handle and illuminator. (HEINE C-145.00.000)

(see page 83).

HEINE® Ophthalmoscopy trainer

Model eye - teaching aid for ophthalmoscopy and use of the (HEINE C-000.33.010)

ophthalmoscope

Refraction scale: -10 to +10 D

Pupillary diameter: variable between 2, 3, 4, 5, 6 and 8 mm

Focal length: F = 18 mm

HEINE® Retinoscopy trainer

Ideal for practising and teaching retinoscopy.

Refraction scale spherical from -7.00 to +6.00 D in steps of 0.5 dpt

Scale: 0° to 180° for cylindrical lenses

Focal length: F = 32 mm

Pupillary diameter: variable between 2, 3, 4, 6 and 8 mm

For re-orders kindly specify the respective no. of the bulb as indicated on the bulb (e.g.: #089)



#57535




Surgical Forceps Titanium



Capsulorhexis Forceps Titanium



Tying Forceps Titanium



Castroviejo Tying Forceps

5 mm shafts, 45 degree angled Length: 108 mm

Girard Tying Forceps

Curved, 4.5 mm tying platform, round handle

Length: 105 mm

Girard Tying Forceps

Straight, tying platform 4.5 mm, round handle

Length: 105 mm

McPherson Tying Forceps

45 degree angled shafts, 4.5 mm tying platform

Length: 115 mm

Tying Forceps

Straight, 4.5 mm tying platform

Length: 85 mm

Tying Forceps

Straight, 4.5 mm tying platform

Length: 115 mm

Stock code: CDQ2-901

Stock code: CDQ2-905


Stock code: CDQ2-906

Stock code: CDQ2-912


Stock code: CDQ2-918


Stock code: CDQ2-919

Different Forceps Titanium



Scissors Titanium





Eye Scissors Titanium

Eye Scissors

Curved, pointed tips, tip to pivot length 24 mm

Stock code: CDQ9-505

Stock code: CDQ9-506

Stock code: CDQ9-601

Stock code: CDQ9-602

Stock code: CDQ9-603

Stock code: CDQ9-604

Stock code: CDQ4-019

Stock code: CDQ10-010

Stock code: CDQ10-012

Stock code: CDQ10-013

Length: 90 mm

Eye Scissors

Curved, pointed tips, tip to pivot length 26 mm

Length: 95 mm

Eye Scissors

Curved, pointed tips, tip to pivot length 32 mm

Length: 110 mm

Eye Scissors

Curved, blunt tips, tip to pivot length 32 mm

Length: 110 mm

Eye Scissors

Straight, pointed tips, tip to pivot length 32 mm

Length: 110 mm

Eye Scissors

Straight, blunt tips, tip to pivot length 32 mm

Length: 110 mm

Specula Titanium

Slade Speculum LASIK

Adjustable, open Length: 65 mm

Liebermann Speculum

Adjustable

V-shaped wires, open

Length: 65 mm

Liebermann Speculum

Adjustable, open, temporal

K-shaped wires Length: 65 mm

Liebermann Speculum

Adjustable, open, nasal

K-shaped wires

it shaped whes

Length: 65 mm



Specula Titanium

Liebermann Speculum

Adjustable, open, nasal

V-shaped wires

Length: 80 mm

Shapira Speculum

Adjustable, open

Length: 65 mm

Speculum

Adjustable, open 7 mm wires

Length: 80 mm

Liebermann Speculum Stock code: CDQ10-011

Stock code: CDQ10-015

Stock code: CDQ10-016

Stock code: CDQ10-030

Stock code: CDQ10-027

Adjustable Solid blades Length: 75 mm

Steinert Speculum Stock code: CDQ10-014

Adjustable Solid blades Length: 80 mm

Genisi Speculum Stock code: CDQ10-031

Adjustable Solid blades Length: 80 mm

Aspiration Speculum

Adjustable Fenestrated Length: 90 mm

Wire Speculum Stock code: CDQ10-020

Length: 52 mm

Barraquer Speculum Stock code: CDQ10-001


Solid blades Length: 52 mm

Needle Holders Titanium



Further instruments on request

Choppers Titanium



Akahoshi Combo Prechopper

Tip with sharp and blunt side Maximum tip opening 3 mm Cross-action, round handle Length: 120 mm Stock code: CDQ2-701

· ·



Ambidextrous, straight Rounded distal tip highly polished Tip 1.1 mm long Length: 115 mm Stock code: CDQ6-515

Drysdale Paddle Manipulator & Phaco Chopper Stock code: CDQ6-520

Double ended

Paddle manipulator, 60 degree phaco chopper

Length: 125 mm

Liebermann Micro Finger & Phaco Chopper Stock code: CDQ6-521

Double ended Length: 125 mm

Phaco Chopper, right Phaco Chopper, left

Unique elongated distal segment, left Special 45 degree inner blade Rounded distal tip highly polished Stock code: CDQ6-524 Stock code: CDQ6-538

Stock code: CDQ6-534

Stock code: CDQ6-536

Stock code: CDQ6-537

Length: 115 mm

Quick Chopper Length: 115 mm

60 degree angled blade, 1.5 mm cutting edge

Length: 115 mm

Nagahara Phaco Chopper Stock code: CDQ6-535

Rosen Phaco Chopper

45 degree angled wedge shaped inferior edge 1.2 mm blunt tip Length: 115 mm

Phaco Chopper, right

Side cutter chopper, right hand Rounded distal tip highly polished Tip 1.8 mm long

Tip 1.8 mm long Length: 115 mm



Manipulators / Hooks Titanium

Synechiae Spatula Stock code: CDQ6-028 Tip 0.5 mm Length: 124 mm Stock code: CDQ6-033 **Drysdale Nucleus Manipulator** Paddle-shaped tip Length: 117 mm **Kuglen Iris Hook & Lens Manipulator** Stock code: CDQ6-048 Straight Length: 115 mm Kuglen Iris Hook & Lens Manipulator Stock code: CDQ6-053 Angled Length: 115 mm Stock code: CDQ6-049 **IOL Manipulator** Angled, curved Length: 126 mm Maltzman-Frenzl Lens Manipulator Stock code: CDQ6-052 Angled V-shaped tip, 0.15 mm diameter Length: 120 mm Sinskey Lens Manipulating Hook Stock code: CDQ6-062 Angled Blunt tip, 0.15 mm diameter Length: 125 mm Jaffe-Bechert Lens Nucleus Rotator Stock code: CDQ6-065 Angled Y-shaped tip, vertical Length: 120 mm **Mushroom Nucleus Manipulator** Stock code: CDQ6-073 Angled Length: 116 mm **Schocket Scleral Depressor** Stock code: CDQ6-070 Double-ended Length: 138 mm

Forceps Stainless steel



Scissors Stainless steel

COO ACT	Vannas Capsulotomy Scissors Straight, tip to pivot length 6.5 mm Length: 90 mm	Stock code: CDQCS9-822
CSO.823	Vannas Capsulotomy Scissors Curved, tip to pivot length 6.5 mm Length: 90 mm	Stock code: CDQCS9-823
	Gills-Vannas Capsulotomy Scissors Straight, tip to pivot length 10 mm Length: 90 mm	Stock code: CDQCS9-824
C39.825	Gills-Vannas Capsulotomy Scissors Curved, tip to pivot length 10 mm Length: 90 mm	Stock code: CDQCS9-825
CO 405	Gills-Vannas Capsulotomy Scissors Angled, tip to pivot length 10 mm Length: 90 mm	Stock code: CDQCS9-826
Company Company	Gills-Vannas Capsulotomy Scissors Angled, curved, tip to pivot length 10 mm Length: 90 mm	Stock code: CDQCS9-827
Car son	Gills-Vannas Capsulotomy Scissors Side-angled, tip to pivot length 10 mm Length: 90 mm	Stock code: CDQCS9-828
3	Vannas Scissors Straight, tip to pivot length 9.5 mm Length: 85 mm	Stock code: CDQCS9-830
Go and	Vannas Scissors Curved, tip to pivot length 9.5 mm Length: 85 mm	Stock code: CDQCS9-831

Eye Scissors Stainless steel



Eye Scissors Stock code: CDQCS9-612

Straight, pointed tips, tip to pivot length 32 mm,

Length: 110 mm

Eye Scissors Stock code: CDQCS9-613

Straight, blunt tips, tip to pivot length 32 mm,

Length: 110 mm

Eye Scissors Stock code: CDQCS9-614

Curved, blunt tips, Length: 105 mm

Eye Scissors Stock code: CDQCS9-615

Straight, pointed tips, Length: 105 mm


Eye Scissors Stock code: CDQCS9-616

Curved, blunt tips, Length: 105 mm

Eye Scissors Stock code: CDQCS9-617

Straight, blunt tips, Length: 105 mm

Specula Stainless steel



Barraquer Speculum for children Stock code: CDQCS10-006

Fenestrated, wires 11 mm, Length: 35 mm

Kratz Barraquer Speculum Stock code: CDQCS10-009

Fenestrated, wires 14 mm, Length: 38 mm

Barraquer Speculum Stock code: CDQCS10-002

Fenestrated, Length: 45 mm

Alfonso Speculum for newborns Stock code: CDQCS10-010

Wires 5 mm Length: 45 mm

Lieberman Speculum Stock code: CDQCS10-021

Adjustable, open, V-shaped wires 15 mm

Length: 85 mm

Kershner Speculum Stock code: CDQCS10-022

Adjustable, open, wires 15 mm

Length: 85 mm

Kershner Speculum for children Stock code: CDQCS10-023

Adjustable, wires 11 mm, Length: 80 mm

Murdock Speculum Stock code: CDQCS10-028

Left side, fenestrated, wires 15 mm

Length: 60 mm

Cook Speculum Stock code: CDQCS10-024

Left side, solid blades 11 mm

Length: 60 mm

Needle Holder Stainless steel




Cannulas Stainless steel











Lacrimal Cannula

Straight, 23G Tip 22 mm long Blunt tip

Irrigating Vectis

Needle 18G Loop 25G Length: 30 mm

Nichamin Hydrodissection Cannula

Angled Curved, 25G Length: 25 mm

Yalon Hydrodissection Cannula

Angled, tip 8 mm Straight, 25G Length: 39 mm

Simcoe Double-Barreled I/A Unit Cannula

Curved, 23G Aspiration port 0.3 mm Length: 28 mm

Simcoe Double-Barreled I/A Unit Cannula Stock code: CDQCS1-035

reserved 23G

Length: 28 mm

Surgical Container trays

Surgical Container tray small 158 mm x 90 mm x 22 mm (L/W/H))

Surgical Container tray medium 225 mm x 142 mm x 22 mm (L/W/H)

Surgical Container tray large 265 mm x 160 mm x 22 mm (L/W/H)

Surgical Container tray double layered 265 mm x 162 mm x 22 mm (L/W/H)

Stock code: CDQ11-009

Stock code: CDQCS1-023

Stock code: CDQCS1-040

Stock code: CDQCS1-026

Stock code: CDQCS1-042


Stock code: CDQCS1-034

Stock code: CDQ11-010

Stock code: CDQ11-001

Stock code: CDQ11-011

SURGICAL INSTRUMENTS



Miscellaneous Titanium

Castroviejo Lacrimal Dilatator Stock code: CDQ7-024

Double-ended

LASIK Flap Elevator Stock code: CDQ4-008

Angled, curved, flattened tip, Length: 70 mm

Mendez Degree Gauge Stock code: CDQ5-205

Round handle, outer Ø 18 mm, inner Ø 12 mm

Mendez Degree Gauge Stock code: CDQ5-211

Round handle, outer Ø 16 mm, inner Ø 12 mm

(without image)

Gravity Marker Stock code: CDQ5-208

For Limbal Relaxing Incisions and toric IOL

implantaion

Graether Half-Circle Axis Marker Stock code: CDQ5-213

Suitable for article 5-205, Mendez Degree

Gauge

Trabeculectomy Punch Stock code: CDQ5-605

0.75 mm

Meyerhoefer Chalazion Curettes

 1 mm, Length: 110 mm
 Stock code: CDQ5-403

 1.5 mm, Length: 110 mm
 Stock code: CDQ5-406

 2 mm, Length: 110 mm
 Stock code: CDQ5-404

 2.5 mm, Length: 110 mm
 Stock code: CDQ5-405

 3 mm, Length: 110 mm
 Stock code: CDQ5-410

 4 mm, Length: 110 mm
 Stock code: CDQ5-411

 4.5 mm, Length: 110 mm
 Stock code: CDQ5-412

Meyerhoefer Chalazion Curettes double ended

 1 mm / 2 mm
 Stock code: CDQ5-402

 2 mm / 4 mm
 Stock code: CDQ5-421

 3 mm / 4 mm
 Stock code: CDQ5-408

Castroviejo Calipers Stock code: CDQ5-053

0 – 20 mm, straight

Braunstein Fixed Calipers Stock code: CDQ5-051

3.5 mm / 4 mm, Length: 70 mm

Halsted Mosquito Forceps Stock code: CDQ2-505

Curved jaws 32 mm, Length: 125 mm

Halsted Mosquito Forceps Stock code: CDQ2-508

Straight jaws 32 mm, Length: 125 mm

Hartmann Mosquito Forceps Stock code: CDQ2-509

Straight jaws 32 mm, Length: 95 mm

Miscellaneous Stainless steel **Wilder Lacrimal Dilatator** Stock code: CDQCS6-104 Length: 115 mm Wilder Lacrimal Dilatator Stock code: CDQCS6-105 Length: 93 mm **Girard Synechiae Spatula** Stock code: CDQCS6-106 Tip 0.5 mm Length: 115 mm Foreign Body Needle Stock code: CDQCS5-006 Straight Length: 120 mm Stock code: CDQCS5-007 **Foreign Body Needle** Curved Length: 120 mm **Phaco Chopper** Stock code: CDQCS6-539 Tip 1.2 mm long Length: 115 mm **Rosen Phaco Chopper** Stock code: CDQCS6-540 60° angled wedge shaped inferior edge, blunt Length: 115 mm **Phaco Chopper** Stock code: CDQCS6-074 Length: 115 mm Frenzl Lens Manipulating Hook Stock code: CDQCS6-076 Length: 115 mm **Kuglen Iris Hook & Lens Manipulator** Stock code: CDQCS6-079 Length: 125 mm **Desmarres Lid Retractors** 6 mm, Length: 138 mm Stock code: CDQCS6-095 8 mm, Length: 138 mm Stock code: CDQCS6-122 11 mm, Length: 138 mm Stock code: CDQCS6-123 13 mm, Length: 138 mm Stock code: CDQCS6-124

15 mm, Length: 138 mm 17 mm, Length: 138 mm Stock code: CDQCS6-125

Stock code: CDQCS6-126

TRU-PACK®

INDIVIDUAL SURGICAL SETS FOR OPHTHALMOLOGY

Surgical sets for immediate use

- for a smooth and efficient procedure
- packed ready to use
- packing sequence precisely tailored to your procedure
- time-saving in emergency surgeries
- standardisation of work processes and reduction of training times



More information at www.tshs.eu/en

Highest quality and best service

TRU-PACK® is a first-class product made in Germany that reflects our know-how and many years of experience in the field of medical devices.



Made in Germany

Since 2016, we have been assembling our customized surgical sets for ophthalmology, which are fitted exclusively with high-quality components and are packaged and sterilized in Germany.

Assembly is carried out by well-trained and experienced employees in our clean room, built especially for this purpose.

ALPHABETIC PRODUCT LISTING

Product	Page	Product	Page
AAPOS Vision - Screening Sets	Page 35	Diagnostic Eye products Miscellaneous	Page 74
Accommodometer	Page 54	Diagnostic Torch Penlight	Page 77
Adhesive Occlusion foils	Page 18	Dinosaur stereo test Bernell	Page 69
Adhesive Prism foils	Page 18	Disposable Occluders	Page 57
Adhesive Spherical foils segments	Page 18	DVD "Leo Learns By Doing"	Page 76
Aesthesiometer Cochet-Bonnet Luneau	Page 74	ETDRS charts	Page 46
Amsler test ring binder	Page 72	ETDRS Illuminated Cabinets	Page 38
Baby trial frames COMO	Page 63	ETDRS Low Contrast charts	Page 50
Bagolini hand-held frame	Page 54	ETDRS Original series – SLOAN letters charts	Page 47
Bagolini lenses	Page 54	Examination Sets	Page 35
Battery handles HEINE®	Page 83	Exophthalmometer	Page 77
Bernell Fixation Stick	Page 22	Eye compresses OPTOMED®	Page 15
Binocularity tests	Page 54 – 61	Eye Models	Page 78
Birkhaeuser reading plates	Page 32	Eye patches ORTOPAD®	Page 6
Block Prism sets	Page 20	Eye Scissors	Page 91
Blue filter for hand-held Slit Lamp	Page 85	Eye Shields	Page 14
Brock strings	Page 77	Film wound dressings	Page 16
C-Test far distance	Page 28	Finoff Transilluminator HEINE®	Page 82
C-Test near vision Haase	Page 35	Fixation cards for HEINE® Retinoscopes	Page 85
Cannulas Stainless steel	Page 100	Fixation Targets	Page 22
Capsulorhexis Forceps	Page 87, 96	Fixed Calipers Braunstein	Page 10
Cardiff Acuity Cards	Page 23	Flashing Fixation Star	Page 77
Castroviejo Calipers Titanium	Page 101	Foldable charts Far distance	Page 28
Chalazion Curettes Titanium	Page 101	Forceps Stainless steel	Page 96
Chalazion Forceps Titanium	Page 89	Forceps Titanium	Page 86
Chargers HEINE®	Page 82	Fundus charts HEINE®	Page 8
Children trial frames	Page 62 – 63	Gravity Marker	Page 10
Choppers Stainless steel	Page 102	Half-Circle Axis Marker Graether	Page 10
Choppers Titanium	Page 94	Hand-held Instruments HEINE®	Page 82
Cilia Forceps	Page 89, 96	HEIDI® Expressions Test Game	Page 25
Clip-on Occluders	Page 57	HEIDI® Fixation Targets	Page 22
Colibri Forceps	Page 89, 96	HEINE® Instruments and accessories	Page 80
Colour tests	Page 70 – 71	HIDING HEIDI® Low Contrast tests	Page 25
COMO Trial Frames children	Page 63	HOTV – foldable chart	Page 37
Confirmation bar	Page 64	HOTV Fixation Cube	Page 22
Confirmation Flippers spherical	Page 65	HOTV Screening flip chart Massachusetts	Page 35
Consumables	Page 73	HRR Colour tests	Page 71
Cover-Test Paddle	Page 56	Illuminated Cabinets and accessories	Page 40
CSV System and charts	Page 38	Indicators for testing VA	Page 78
Damato Campimeter	Page 72	Individual Prisms	Page 19
Desmarres Lid Retractors	Page 102	Injection plasters	Page 16

ALPHABETIC PRODUCT LISTING

Product	Page	Product	Page
Ishihara Colour tests	Page 70	Maddox Cross	Page 61
Jackson Cross cylinders	Page 64	Maddox Cylinder for use in trial frames	Page 54
Jaeger reading card	Page 32	Maddox Cylinders	Page 56, 66
KAY PICTURES	Page 36	Maddox-Wing test	Page 61
Lacrimal Dilatator Titanium	Page 101, 102	Manipulators/Hooks Stainless steel	Page 95, 102
Lambda 100 Retinometer HEINE®	Page 82	MARS Low Contrast tests	Page 50
Landolt-C Response items	Page 78	Marsden Fixation ball	Page 74
Landolt-C charts translucent	Page 41	Mendez Degree Gauge	Page 101
Landolt-C foldable chart	Page 28	Metal Rim lenses Overview	Page 66
LANG Fixation Stick and Cube	Page 22	Mini 3000® (LED Ophthalmoscope HEINE®)	Page 80
LANG Stereo tests	Page 68	Mirror Exophthalmometer Hertel	Page 77
Lanthony colour test (desaturated Farnsworth D-15)	Page 70	Mosquito Forceps	Page 101
LASIK Flap Elevator	Page 101	Multiple Pinhole Spectacle Occluder	Page 57, 75
LEA 3-D PUZZLE®	Page 25	Near Low Contrast tests – ring binder	Page 53
LEA™ Book "What and how does this child see?"	Page 76	Near vision tests	Page 32 – 35
LEA [™] Cards Low Contrast ring binder	Page 53	Needle Holders	Page 93, 99
LEA™ Cone Adaptation test	Page 24	Neuro-Orthoptic tests	Page 26 – 27
LEA™ Core Tests - set	Page 29	Nose pads for trial frames	Page 75
LEA™ Domino Set	Page 24	Nystagmus glasses Frenzel	Page 61
LEA™ Fixation Targets	Page 22	Occluders and Paddles	Page 56 – 57
LEA GRATINGS® Acuity Test and Paddles	Page 26	Occlusion glasses	Page 58
LEA GRATINGS® Contrast Test	Page 22, 52	OKN Instruments	Page 60
LEA™ Low Contrast charts	Page 45	Opaque Glass Occluders for children	Page 58
LEA MAILBOX GAME®	Page 24	Ophthalmic Examination Lamp HEINE®	Page 82
LEA™ Panel 16 – Colour tests	Page 71	Ophthalmology Sets HEINE®	Page 81
LEA RECTANGLES GAME®	Page 24	Ophthalmoscopes HEINE®	Page 80
LEA™ Ring Binders	Page 29, 53	Ophthalmoscopy Trainer HEINE®	Page 85
LEA™ Ruler Pads	Page 78	OPTOMED® Eye compresses	Page 15
LEA™ Single Presentation Cards	Page 25	Orange filter for BETA 200® retinoscope	Page 85
LEA™ Stereo tests	Page 69	ORTOLUX®/ORTOLUX® AIR eye shields	Page 14
LEA SYMBOLS®/ Tumbling E-hooks Wallchart	Page 28	ORTOPAD® Eye patches	Page 6 – 10
LEA SYMBOLS®/Landolt rings Wallchart	Page 28	Paediatric Vision Tests	Page 21
LEA SYMBOLS® Playing Cards	Page 24	Paediatric Visual Acuity Set	Page 35, 58
LEA SYMBOLS® – Screening Sets	Page 33	Patient card LAMBDA 100 Retinometer HEINE®	Page 82, 85
LEA™ Translucent charts	Page 48	Polarisation Glasses	Page 75
LEA™ USB stick "What and How Does This Child See?	" Page 76	Polarisation Occluders for children	Page 58
Lid retractors Desmarres	Page 102	Polarisation/Red-Green Flippers	Page 57
Loupe attachment for hand-held Slit Lamp HEINE®	Page 85	Preferential Looking tests	Page 23
Low Contrast Test charts and ring binders	Page 50 – 53	Prism Bars, Prism Sets and accessories	Page 19, 20
Low Vision tests	Page 30 – 31	Prism Flippers	Page 65
LUVIS mini – Illuminated Cabinet	Page 39	Provikit eye tests (Project Provision)	Page 35, 58

ALPHABETIC PRODUCT LISTING

		_	
roduct		Page	
Pupillometer		Page 77	
Reading charts – hiç	gh/10 % low contrast	Page 34	
Reading plates Low V	/ision	Page 53	
Rechargeable Batte	ries HEINE®	Page 83	
Recording forms for	colour tests	Page 70	
Red-Green glasses,	frames & flippers	Page 55	
Red-Green tests		Page 54	
Refraction and acce	ssories	Page 62 – 67	
Response plates an	d cards	Page 76	
Retina Lenses Medi\	Norks	Page 67	
Retinometer HEINE®	LAMBDA 100	Page 82	
Retinoscopes HEINE	BETA® 200	Page 81	
Retinoscopy Bar sets	5	Page 64	
Retinoscopy trainer	HEINE®	Page 85	
Reversible Occluding	g Glasses	Page 76	
RUNGE near vision of	cards	Page 32	
Scissors		Page 90, 97	
Single-use Tonomet	er Prisms L900 Luneau	Page 73	
SLOAN letters – folde		Page 28	
SLOAN Low Contrast	charts	Page 45	
Spare bulbs HEINE®		Page 84	
Spectacle Occluders	;	Page 57, 75	
Specula Stainless steel		Page 98	
Specula Titanium		Page 91 – 92	
ipielmann Occluder	•	Page 56	
Spool Plasters		Page 15	
Stereo tests		Page 68 – 69	
Surgical Container T	rays	Page 100	
urgical Instruments		Page 86 – 102	
Surgical Instruments	s Miscellaneous	Page 101 – 102	
Feller Acuity Cards®	II	Page 23	
Translucent charts D	istance, near, intermediary distance	Page 41 – 49	
Translucent Low Cor	•	Page 45	
Trial Frames		Page 62 – 63	
Trial Lenses set		Page 64	
Tumbling E Respons	e optotype	Page 78	
Tying Forceps		Page 88, 96	
Visual Development	/Rehab tests	Page 24 – 25	
Visual Field Fatigue		Page 57, 73	
Visual Field tests		Page 72	
	eSt "Color Vision Testing Made Easy"	Page 71	



TRUSETAL VERBANDSTOFFWERK GMBH

Our brand EYESFIRST® offers you a comprehensive range of GOOD-LITE and LEATM test products commonly used in Europe as well as products of other highly renowned manufacturers.

In addition to the catalogue, we kindly request you to visit our online shop **www.eyesfirst.eu**, which will be continuously updated during the course of the year.

We are happy to answer your requests any time.



Rita Brüggemann CEO



Jessie Ilsen CEO



Tanja Kerins Head of TRU-PACK®



Ute Rau
Head of EYESFIRST®



The company office of Trusetal in Schloss Holte-Stukenbrock (Germany). Your specialist for all your Ophthalmic & Orthoptic products including surgical instruments, dressing materials and wound care needs.



We look forward to being of service to you! Should you have any questions, suggestions or comments please contact us by phone, fax or email.

Looking forward to your call and your visit to our website!

OUR RENOWNED BRANDS:

EYESFIRST®
ORTOPAD®
TRU-PACK®

TRUSETAL VERBANDSTOFFWERK GMBH

Konrad-Zuse-Str. 15 33758 Schloss Holte-Stukenbrock Germany

Phone: +49 5207 991688-0

Fax: +49 5207 991688-28

E-Mail: info@eyesfirst.eu

info@tshs.eu

Internet: www.tshs.eu/en

www.augenpflaster.com

Our online shop:

www.eyesfirst.eu

Our email address for your orders:

orders@eyesfirst.eu